NEW BUCHI Kjeldahl Tablets





The Kjeldahl Tablets

- Massive acceleration of Kjeldahl digestion
- Consist of SULFATE and METAL SALTS



BUCHI Range of Kjeldahl Tablets

	Article (box of 1000 pcs.)	Benefit	Composition
	Titanium	Time saving	3.5 g K ₂ SO ₄ / 0.105 g CuSO ₄ · 5 H ₂ O
	# 11057980		0.105 g TiO ₂
	Titanium Micro	Reduced	1.5 g K ₂ SO ₄ / 0.045 g CuSO ₄ · 5 H ₂ O
	# 11057981	chemical amount	0.045 g TiO ₂
	Missouri	Easy to use and	4.98 g K ₂ SO ₄
	# 11057982	universally applicable	0.02 g CuSO ₄ • 5 H ₂ O
	ECO	Environmentally	3.998 g K ₂ SO ₄
	# 11057983	friendly	0.002 g CuSO ₄
	Antifoam	Maximum foam	0.97 g Na ₂ SO ₄
	# 11057984	reduction	0.03 g Silicone Antifoam
	Copper Micro	Reduced	1.5 g K ₂ SO ₄
	# 11057985	chemical amount	0.15 g CuSO ₄ • 5 H ₂ O

Choosing the Correct Tablet

What to consider when choosing a tablet?

.....

Ideal digestion conditions
are achived with
2 ml H₂SO₄ to 1 g of catalyst

- Boiling point of 370°C
- No loss of nitrogen
- Minimal time needs

Recommendation

Weight

- Safety aspects
- Temporal aspects
- Ecological aspects
- Foam formation
- Chemical amount

Use the enclosed decision wheel to select the tablet that best suits your needs.

3.71 g	Optimal compromise between environmental and performance priorities.
1.59 g	Same as Titanium (11057980) but for semi-micro-Kjeldahl & micro-Kjeldahl applications.
5 g	The digestion with the Missouri catalyst is more environmentally friendly.
4 g	Our most environmentally friendly catalyst, due to the very low copper content.
1 g	Used as a general purpose foam suppressant. This catalyst has to be combined with Titanium Micro (11057981) or Copper Micro (11057985).
1.65 g	Combo tablets for antifoam or for micro Kjeldahl applications

Use the decision wheel to choose the catalyst that best suits your needs.


- 1. Select the correct side of the wheel according to your application
- Standard Kjeldahl method
- Micro-Kjeldahl method
- 2. Rotate the wheel till your decision factors appear
- Foam formation of sample
- Time requirements for digestion
- Environmental importance

Example

Fast digestion of a non-foaming sample

- 1. Standard Kjeldahl application
 - select side for "standard Kjeldahl"
- 2. No foam formation
 - 4 foam formation "no"
- 3. Fast digestion
 - 4 temporal aspects important "2"

Recommended tablet: Titanium



BÜCHI Labortechnik AG CH-9230 Flawil 1

T +41 71 394 63 63

F +41 71 394 65 65

buchi@buchi.com

www.buchi.com/kjeldahl-tabs

