MultivaporTM Application Guide





en

Table of Contents

This document describes all aspects which are important for the evaporation of various solvents using the Multivapor. This handbook is a guideline for the evaporation of your specific solvent and mixtures at high performance. Some standard applications are described in detailed Application Notes. best@buchi are found on the internet under www.buchi.com, then Applications and best@buchi. Please contact your local BUCHI representative for additional information on a particular application.

| 1 | Introd | luction | 3 |
|---|--------|---|----|
| | 1.1 | Multivapor and Configurations | 3 |
| | 1.1.1 | Standard configuration | 3 |
| | 1.1.2 | Variations and Versions | 4 |
| | 1.1.3 | Recommended accessories | 4 |
| | 1.1.4 | Limits of the Multivapor | 4 |
| | 1.2 | The process of distillation | 5 |
| | 1.3 | General guidelines for evaporation of solvents | 5 |
| | 1.3.1 | Correct load of the condenser | 6 |
| | 1.3.2 | Manual distillation and the Δt 25/20°C rule | 6 |
| | 1.3.3 | Distillation with gradients | 6 |
| | 1.3.4 | EasyVac function | 7 |
| | 1.4 | Parameters and impact on distillation | 8 |
| | 1.4.1 | Heat medium | 8 |
| | 1.4.2 | Vacuum | 8 |
| | 1.4.3 | Temperature of condenser | 8 |
| | 1.4.4 | Rotation speed | 9 |
| | 1.4.5 | Vessel type | 9 |
| | 1.4.6 | Solvent volume | 9 |
| | 1.4.7 | Solvent mixtures | 9 |
| | 1.4.8 | Multivapor P-6 and P-12 | 10 |
| | 1.4.9 | Troubleshooting | 10 |
| | 1.5 | Sample applications | 11 |
| 2 | Evapo | pration of pure solvents | 13 |
| | 2.1 | List of solvents and classification | 13 |
| | 2.2 | Calculation of boiling point from applied pressure | 13 |
| | 2.3 | Evaporation of low-boiling solvents | 14 |
| | 2.3.1 | Solvent: Acetone (≥99.8%) | 14 |
| | 2.3.2 | Solvent: Ethyl acetate (Analytical grade) | 15 |
| | 2.3.3 | Solvent: <i>n</i> -Hexane (purum ≥95.0%) | 16 |

| | 2.3.4 | Solvent: <i>n</i> -Hexane – P-6 (purum ≥95.0%) | 17 |
|---|-----------------------------------|---|----------------------------|
| | 2.4 | Evaporation of medium-boiling solvents | 18 |
| | 2.4.1 | Solvent: Xylene (Mixture of isomers, purum ≥97.0%) | 18 |
| | 2.5 | Evaporation of high-boiling solvents | 19 |
| | 2.5.1 | Solvent: Dimethylformamide, DMF (purum, ≥98.0%) | 19 |
| | 2.6 | Evaporation of polar protic solvents | 20 |
| | 2.6.1 | Solvent: Ethanol (4.8% methanol, purum, ≥99.8%) | 20 |
| | 2.6.2 | Solvent: Distilled water, H ₂ O | 21 |
| 3 | Evapo | pration of solvent mixtures | 22 |
| | | | |
| | 3.1 | Hexane / dichloromethane 7:3 (e.g. flash chromatography) | 22 |
| | 3.1 3.2 | Hexane / dichloromethane 7:3 (e.g. flash chromatography) Ethanol / water 6:4 (e.g. reverse phase-chromatography) | 22 24 |
| | | | |
| 4 | 3.2 3.3 | Ethanol / water 6:4 (e.g. reverse phase-chromatography) | 24 |
| 4 | 3.2 3.3 | Ethanol / water 6:4 (e.g. reverse phase-chromatography) Azeotrope chloroform / methanol 8:1 (e.g. flash chromatography) | 24 26 |
| 4 | 3.2 3.3 Apper | Ethanol / water 6:4 (e.g. reverse phase-chromatography) Azeotrope chloroform / methanol 8:1 (e.g. flash chromatography) ndix A – Reference BUCHI Application Notes | 24 26 29 |
| 4 | 3.2 3.3 Apper 4.1 4.2 | Ethanol / water 6:4 (e.g. reverse phase-chromatography) Azeotrope chloroform / methanol 8:1 (e.g. flash chromatography) ndix A – Reference BUCHI Application Notes Fat Determination in Shortbread using SpeedExtractor E-916 | 24 26 29 29 |
| - | 3.2 3.3 Apper 4.1 4.2 | Ethanol / water 6:4 (e.g. reverse phase-chromatography) Azeotrope chloroform / methanol 8:1 (e.g. flash chromatography) ndix A – Reference BUCHI Application Notes Fat Determination in Shortbread using SpeedExtractor E-916 Flash chromatography on cartridges for the separation of plant extracts | 24 26 29 29 31 |

© BÜCHI Labortechnik AG, 2012, Dr. D. Rütti, M. Heierli

1 Introduction

1.1 Multivapor and Configurations

The Multivapor is a compact 6 or 12 position parallel vortex evaporator. Since the efficiency of evaporation depends on a variety of parameters, careful optimization for each type of sample/evaporation task is important. This handbook provides practical guidelines to successfully develop new applications with the goal that you could enjoy your Multivapor.

1.1.1 Standard configuration

If not indicated otherwise, the following Multivapor configuration was used in this booklet (**Figure 1**).

| Component | Configuration | | |
|--|--|--|--|
| Multivapor | P-12 (MP11001S12) 12 positions | | |
| Condenser | PETP Cover and Adapters Type S condenser for tap water or a recirculating chiller | | |
| Vacuum controller V-855, Software V 3.03 | | | |
| Vacuum pump | V-700, Woulff bottle, Secondary condenser | | |
| Chiller | F-105, Set to 5°C Cooling power at 15°C: 500 W Flow rate (no accessories): 1.8 l/min | | |
| Vessel | BUCHI SpeedExtractor tube 60 ml 30 ml max working volume (49535) | | |
| Options | None (no Woulff bottle, no Protective shield, no frits) | | |



Figure 1: BUCHI Multivapor P-12 with type S condenser, vacuum pump V-700 and vacuum controller V-855

1.1.2 Variations and Versions

| Component | Variants |
|-------------------------------------|------------------------------------|
| Tube adapter | PETP (white, blind adapter: black) |
| | PEEK (ecru) |
| Seals for tube adapters | PTFE |
| Vacuum cover | PETP (black) |
| | PEEK (ecru) |
| O-ring for the vacuum cover adapter | EPDM (black) |
| | FKM (green) |
| | FFKM (black) |
| Vacuum cover seal | EPDM |
| Vacuum tube | Ribbed PFA with PTFE seals |

The PETP vacuum cover is equipped with EPDM O-rings. A set of 12 FKM O-rings is enclosed. Optionally, there are highly resistant FFKM O-rings available. The PEEK vacuum cover is equipped with FFKM O-rings by default.

For order codes please consider the operation manual or contact your local dealer.

1.1.3 Recommended accessories

The most convenient way to achieve the right setting is to use the BUCHI vacuum controller V-855 (**Figure 2**). With this instrument, the pressure is either determined automatically using the EasyVac function (P-12 only), or it can be easily determined by means of an integrated temperature/pressure converter (the so-called Solvent Library) which is available for a variety of solvents. Moreover, pressure gradients can be programmed to ensure smooth and reproducible evaporation with the highest possible yields.¹

The use of a chiller (**Figure 3**) is highly recommended. A chiller enables a constant low temperature without tap water input. A low condenser temperature is needed to distill at low temperatures without loss of solvents into the environment. BUCHI chillers can be controlled by the V-855 vacuum controller.

1.1.4 Limits of the Multivapor

The limit is given by the lowest sensible working pressure of approximately 25-30 mbar (V-700) and the unheated cover leading to condensation in the cover of certain solvents. For high-boiling solvents, the Syncore Polyvap with heated cover is recommended.



Figure 2: BUCHI vacuum controller V-855 (professional)



Figure 3: BUCHI recirculation chiller

¹ Optimizing the evaporation using the Syncore® Analyst, recovery of 94±3% was achieved for naphthalene in PAH determination of soil samples. See: J. Cavegn, B. Haag, R. Hartmann, Trace analysis of PAHs: evaluation of two extraction methods (EPA3541 and EPA3545) and optimization of subsequent concentration, best@buchi no. 51/2008.

The Multivapor is limited when the task is to concentrate to a defined volume instead to dryness. For this application, the Syncore Analyst with cooled appendix was developed.

1.2 The process of distillation

Vacuum evaporation is the process of causing the pressure in a liquidfilled container to be reduced, causing the liquid to evaporate at a lower temperature than normal. A distillation consists of vaporization and subsequent condensation. When these two processes are well parted, separation of substances is possible based on differences in vapor pressure. The transport of the vapor is driven by a thermal pump, a dynamic difference in pressure between places of evaporation and condensation.

The evaporation speed is positively correlated with the amount of heat supplied by the heating plate. However, there is no correlation between the set temperature on the instrument (T_i) and the temperature of the sample (T_s). Liquids cannot assume a temperature higher than their boiling point (T_{bp}), which is solely determined by the applied vacuum. The only critical point for the sample is at the end of an evaporation to dryness process (**Figure 4**).

In the condenser, the heat of evaporation is then withdrawn from the vapor; the vapor condenses.² For best performance and highest recovery of the solvent, it is important to maintain a balanced evaporation-condensation process. Whenever more substance is vaporized than is condensed in the cooler, the pressure within the apparatus rises. To maintain the set pressure, the vacuum pump draws continually pulling the vaporized solvent into the environment. This must be avoided.³ Pollution is prevented by assembling a secondary condenser behind the pump, or preferred, adjusting the evaporation to the cooling capacity of the condenser. Ideal conditions are achieved when the condensate covers approximately the half of the condenser (**Figure 5**).

1.3 General guidelines for evaporation of solvents

To achieve the best evaporation conditions, five parameters, the temperature of the heating plate, the boiling point, the cooling water temperature, and the rotational speed of the orbital movement, must be attuned to one another. A tight system enables best distillation performance. Therefore, the apparatus should be checked for leaks without the solvent in it. A well-sealed system keeps the vacuum losses within bounds (< 10 mbar / 5 min at vacuum of 20 mbar).




Figure 4: Change of the sample temperature during the evaporation process with a sudden rise in temperature at the end of the process.

² This phenomenon is referred to as a thermal pump. See: The laboratory assistant, BÜCHI Labortechnik AG, second edition, p. 24.

³ It is estimated that worldwide approx. 160'000 liters of solvent per day are discharged into the air due to unbalanced evaporation and condensation speeds. See: The laboratory assistant, BÜCHI Labortechnik AG, second edition, p. 47.

1.3.1 Correct load of the condenser

Evaporated solvent is condensed in the condenser and evaporation heat is dissipated. Best performance and highest recovery of the solvent are generally achieved when the condensate covers approximately the half of the height of the condenser (**Figure 5**).

1.3.2 Manual distillation and the Δt 25/20°C rule

A simple rule of thumb can be used as a starting point to find good settings for the heating medium, the vapor temperature, and the cooling temperature: the Δt 25/20°C rule. This rule specifies the temperature difference between the three zones. For instance, when the temperature on the instrument is set to 55°C, the vacuum should be set that the boiling point is 30°C, which is 25°C lower. In order to achieve sufficient condensing, cooling temperature of another 20°C lower also 10°C is recommended (**Figure 6**).

The following steps are performed in practice:

- 1. Set the heating medium to the desired temperature according to the stability of the product and the solvent
- Determine the pressure using the Solvent Library so that a boiling point of 25°C lower results
- 3. Set the cooling temperature another 20°C lower

If the distillation performs too slowly, either increase vacuum or heating temperature. Aerate respectively decrease vacuum when sample is foaming or bumping. Distillation of problematic solutions can be distilled without supervision by using the gradient function (see Section 1.3.3). Further troubleshooting is found in detail in Section 1.4.9.

1.3.3 Distillation with gradients

Distillation of problematic solutions can be handled without supervision by using the gradient function. Setting the pressure manually and immediately to the corresponding boiling point accommodates the risk of bumping. To avoid boiling retardation or loss of analyte, program a pressure gradient. Gradients also allow convenient integration of "drying" steps by applying a very low pressure during some minutes after the effective evaporation. To prevent reevaporation of solvent from the receiving flask during "drying", cool the distillate with ice bath or refrigerated receiver.

Optimization of the conditions for your specific problem is usually started with manual distillation. Define the desired boiling point according to the sensitivity of the product and the Δt 25/20°C rule. Start 500 mbar above the calculated pressure. Decrease the pressure quickly by 300-400 mbar, and then approximate slowly the desired




Figure 5: The vacuum and cooling conditions are ideal when the condensate covers approximately half of the height of the condenser.




Figure 6: Illustration of the Δt 25/20°C rule.

boiling point. Keep the pressure constant until the evaporation is finished. For highly volatile compounds slow aeration over 1-2 min is recommended.

| Solvent | | <i>n</i> -hexane | 800 1 |
|----------------------|----------------------|------------------|---------------------------------|
| Desired boiling poin | it, T_{bp} | 30°C | 700 006 mpar 000 |
| Pressure derived fro | om the solvent list | 241 mbar | |
| Instrument tempera | ture, T_i | 55°C | |
| Pressure gradient | | | e 400 300 200 100 0 |
| Starting point | 241 + 500 = 741 mbar | | 0 |
| Ramp | 741-291 mbar | 2 min | |
| | 291-237 mbar | 2 min | |
| | 237-230 mbar | 2 min | |
| | 230-202 mbar | 3 min | |
| Constant | 202 mbar | 1 min | _ |





A first rapid decrease in pressure prevents heating up the sample. The following gradient decreases bumping, retardation and foaming of the product. After a satisfying manual distillation, the repeat function (V-855 only) displays the last distillation for easy documentation and repetition.

1.3.4 EasyVac function

The EasyVac mode (Vacuum controller V-855, software version 3.03 or higher) serves to carry out an automatic and gentle distillation based on the vapor pressure of the solvent. The process measures pressure differences over time and adjusts the vacuum automatically depending on the efficiency of the condenser. This mode is suited for distilling single solvents with the Multivapor under certain prerequirements.

| Prerequirement | | | | |
|-----------------------------|---------|--|--|--|
| Rack temperature | ≥ 50°C | | | |
| Solvents with boiling point | < 100°C | | | |
| Total evaporated volume | ≥ 80 ml | | | |
| Only weak foaming products | | | | |
| Glassware diameter | < 12 mm | | | |

The operator may interact at any time during automatic distillation and the process continues according to the manual set parameters.

| EasyVac Multivapor P-12 | Detection start point | Distillation | Detection end point |
|-----------------------------------|--------------------------|--------------|------------------------|
| Glassware diameter 13 < x < 19 mm | \checkmark | \square | |
| Glassware diameter > 19 mm | \checkmark | \checkmark | V |

| EasyVac Multivapor P-6 | Detection start point | Distillation | Detection end point |
|------------------------|--------------------------|--------------|------------------------|
| Low-boiling solvents | V | (図) | Often too early |

Please note: A tight system is crucial for correct determination. The EasyVac mode only works with a speed controlled vacuum pump (V-700 or V-710). The distillation time can be increased by 30 to 100% compared to manual distillation.

1.4 Parameters and impact on distillation

1.4.1 Heat medium

Supplying more energy to the solvent by elevated heat medium temperature will accelerate the distillation process. This additional energy will dissipate in the condenser leading to possible losses of solvent to the environment due to an overburdened condenser. Leaking solvent might condense in the pump leading to decreased pump performance. Acceleration of the distillation increases risk of foam and bumping sample. Another danger is the loss of analytes due to increased temperature at the end of an evaporation process.

Water in the crystal rack ensures good heat transfer to the sample. Best evaporation performance is achieved when the water level reaches at least the solvent level. However, do not overfill or spill of water by the vortex is expected.

1.4.2 Vacuum

Similarly, decreasing of boiling point of the solvent by reduction of pressure will accelerate the distillation process. This additional energy will dissipate in the condenser leading to possible loss of solvent in the environment due to an overburden of the condenser capacity. If the system is insufficiently tight for the applied vacuum, the pump continually draws solvent, leading to known problems. Analytes might be lost due to boiling. Acceleration of the distillation increases risk of foam and bumping sample.

1.4.3 Temperature of condenser

For efficient condensation of the vapor, the temperature of the condenser should be at least 15-20°C less of the vapor temperature. If the temperature of the cooling medium is fixed e.g. by the temperature of tap water, adjust the pressure and vapor temperature accordingly higher. Otherwise, exhaust of vapor into the environment might occur. Ideal conditions are achieved when the condensate covers approximately the half of the condenser.

Increase distillation speed Increase temperature of heat

medium **Risks** Loss of analytes Heat degradation of sample Foam and bumping of sample

Increase distillation speed Increase vacuum Risks Loss of analytes Foam and bumping of sample

1.4.4 Rotation speed

The surface area of the sample influences the evaporation. The larger the surface, the faster is the evaporation. The surface area is determined by the rotational speed of the heating plate. A faster rotation increases the surface area and hence accelerates the process. However, increasing the rotational speed also increases the vortex of the heating medium water. This might result in a reduced heat transfer surface.

Generally speaking, start the optimization with a pronounced vortex. For the Multivapor, this is typically the case at position 8. Then increase rotational speed until an optimal position for evaporation speed and handling of the sample (foam, bumping) is found.

| 1.4.5 Vessel ty | /pe |
|-----------------|-----|
|-----------------|-----|

The vessel shape (**Figure 7**) influences greatly the distillation performance. Especially the length and diameter are important characteristics. Distillation speed of ethanol and acetone were decreased by approx. 50% in a 50 ml Falcon tube (29 x 114 mm, conical end) and reaction tube (16 x 130 mm) compared to the SpeedExtractor tube (28 x 140 mm). Applied pressures could be transferred but the gradient time will be extended accordingly.

The vessel for quickest distillation is as long and as wide in diameter as possible without conical parts.

1.4.6 Solvent volume

The volume of evaporated ethanol showed no effect on the speed using SpeedExtractor tubes. However, highly volatile solvents and large volumes might lead to starvation of the heat medium, thus elongating the evaporation. A further increase of the vacuum is advisable.

P-12 compatible tubes are preferably filled to only 50% of the total volume with sample. Standard P-6 tubes and SpeedExtractor vessels have working volumes of 150 and 180 ml, respectively.

Applied pressures could be transferred but the gradient time will be extended accordingly.

1.4.7 Solvent mixtures

Mixtures are evaporated by removing one solvent after the other in sequence according to the increasing boiling points. The purpose of this, in general, is not to achieve a sharp separation of the individual components, but rather to concentrate on the dry residue. A number of problems such as boiling delays or undesirable solvent emissions can occur in the area of transitioning from one component to the next.⁴ In addition to this, a time or concentration dependent foam build-up can occur.

| In | crease distillation speed | | | |
|---------------------|----------------------------|--|--|--|
| Rotation speed 8-10 | | | | |
| R | isks | | | |
| S | pill of heat medium | | | |
| L | arge vortex decreases heat | | | |
| tr | ansfer | | | |
| F | oam and bumping of sample | | | |



Figure 7: Selection of compatible P-12 tubes. *From left*. BUCHI SpeedExtractor tube, 50 ml Falcon and test tube.

⁴ A. Talamona, G. Schönenberger, Vacuum control during vaporization of problematic solutions, best@buchi no. 06/2001.

These dangers are avoided if the system pressure is reduced slowly using the gradient function. Start with optimal conditions for lowerboiling solvent according to the gradient rule.

Some mixtures of a certain composition, however, cannot be separated by distillation. They are called azeotrope. When an azeotrope is boiled, the resulting vapor has the same ratio of constituents as the original mixture.

Do not go under lowest condensation temp with boiling point of lowestboiling solvent, or use a well-cooled receiving vessel. Otherwise, reevaporation from receiving flask might occur reducing evaporation performance.

1.4.8 Multivapor P-6 and P-12

Using SpeedExtractor vessels for P-6 and P-12, distillation speed increased by 20% with P-6 compared to P-12 (tested with *n*-hexane). Applied pressures could be transferred but the gradient time will be extended accordingly.

1.4.9 Troubleshooting

| Problem | Problem Action | | |
|----------------------------|--|---|--|
| Foam | Reduce vacuum Program ramp Install frit in the adapter | Section 1.4.2 Section 1.3.3 See Multivapor manual | |
| Overheating, Bumping up | Reduce vacuum Program ramp Reduce heat Use vessel with larger diameter Reduce time span from immersing the vessels into the heat medium to start vacuum and rotation | Section 1.4.2 Section 1.3.3 Section 1.4.1 Section 1.4.5 | |
| Condensation in cover | Reduce boiling temperature (increase vacuum) and adjust heat medium as well as condenser temperature accordingly Reduce length of vacuum tubing to condenser | | |
| Slow distillation speed | Improve heat transfer Optimize rotational speed Increase vacuum Increase temperature | Sections 1.4.1, 1.4.4, 1.4.8, 1.4.5 Section 1.4.4 Section 1.4.2 Section 1.4.1 | |

1.5 Sample applications

The Multivapor was designed to fit excellently in your analytical work flow without transferring of vessels. Samples are often concentrated after extraction (BUCHI Extraction units, SpeedExtractor) or reaction processes as well as after separation by flash chromatography (BUCHI Sepacore).

Extraction: Fat Determination in Shortbread

The determination of fat in food and feed is a routine procedure for quality assurance and labeling. Fat was determined in butter shortbread after homogenization of the sample and extraction with the SpeedExtractor E-916. The solvent was evaporated in parallel using the Multivapor P-6. The total fat content was determined gravimetrically after the extract has been dried to a constant weight.⁵



Figure 8: Shortbread

Chromatography: Separation of plant extracts

Flash chromatography systems can be used not only for rapid and easy purification of synthetic products, but also for the separation of complex natural product mixtures such as plant extracts.

Sepacore cartridges enabled a good separation of compounds with a broad range of polarity, as typically found in plant extracts. The chromatographic resolution remained, however, lower than that achieved by medium pressure liquid chromatography (MPLC) on columns packed with material of smaller particle size. For poorly soluble extracts, solid introduction gave better results than liquid injection.⁶

To recover these extracts in pure form, the mobile phase of the selected fractions will be evaporated. This is preferrably performed using Multivapor, as the collection vessel can be used directly for evaporation without sample transfer.



Figure 9: Separation of plant extracts

⁵ Fat in Shortbread using SpeedExtractor E-916 and Multivapor P-6. BUCHI Application Note no. 005/2009.

⁶ Weber P, Schafroth N, Hamburger M, Potterat O. 2010. Flash chromatography on cartridges for the separation of plant extracts – Rules for the selection of chromatographic conditions, and comparison with MPLC. Poster in collaboration with University of Basel. http://www.buchi.com/Applications.19910.0.html.

2 Evaporation of pure solvents

2.1 List of solvents and classification

| Solvent | Enthalpy of Boiling vaporization point | | Density [g/ml] | Classification | Constant b | |
|--------------------------------|--|--------|-------------------|----------------|----------------|-------|
| | • | [J/ml] | [J/g] | | | |
| Acetic acid | 118°C | 729 | 695 | 1.049 | Polar protic | 0.183 |
| Acetone | 56°C | 437 | 553 | 0.79 | Low boiling | 0.196 |
| Acetonitrile (MeCN) | 82°C | 570 | 725 | 0.786 | Low boiling | 0.195 |
| Benzene | 80°C | 481 | 548 | 0.877 | Low boiling | 0.202 |
| n-Butanol | 118°C | 502 | 620 | 0.81 | Polar protic | 0.155 |
| Chlorbenzene | 132°C | 417 | 377 | 1.106 | Medium boiling | 0.202 |
| Chloroform | 61°C | 392 | 264 | 1.483 | Low boiling | 0.203 |
| Cyclohexane | 81°C | 303 | 389 | 0.779 | Low boiling | 0.206 |
| Cyclopentane | 40°C | 313 | 417 | 0.751 | Low boiling | 0.207 |
| Dichloromethane (DCM) | 40°C | 439 | 330 | 1.33 | Low boiling | 0.194 |
| Diethyl ether | 35°C | 278 | 389 | 0.714 | Low boiling | 0.200 |
| Dimethylformamide (DMF)* | 153°C | 549 | 578 | 0.949 | High boiling | 0.180 |
| 1,4-Dioxane | 101°C | 400 | 388 | 1.033 | Medium boiling | 0.200 |
| Ethanol | 78°C | 694 | 879 | 0.789 | Polar protic | 0.159 |
| Ethyl acetate | 77°C | 355 | 394 | 0.9 | Low boiling | 0.189 |
| Formic acid | 101°C | 601 | 493 | 1.22 | Polar protic | 0.200 |
| n-Hexane | 69°C | 243 | 368 | 0.66 | Low boiling | 0.206 |
| Isopropanol (IPA) | 82°C | 549 | 699 | 0.786 | Polar protic | 0.154 |
| Methanol | 65°C | 971 | 1227 | 0.791 | Polar protic | 0.167 |
| Pentane | 36°C | 239 | 381 | 0.626 | Low boiling | 0.214 |
| n-Propanol | 97°C | 633 | 787 | 0.804 | Polar protic | 0.154 |
| Tetrahydrofuran (THF) | 66°C | 395 | 444 | 0.889 | Low boiling | 0.192 |
| Toluene | 111°C | 370 | 427 | 0.867 | Medium boiling | 0.202 |
| Water | 100°C | 2261 | 2266 | 0.998 | Polar protic | 0.167 |
| Xylene (mixture of isomers) | 138.5°C | 338 | 389 | 0.87 | Medium boiling | 0.199 |

*) Distillation of high-boiling solvents is not recommended with the Multivapor.

2.2 Calculation of boiling point from applied pressure


$$T_p = \frac{T_{bp}}{(3.006 - \log p) \cdot b + 1}$$
 (Equation 1)

 T_p [K] is the temperature of the boiling point at the pressure *p* [mbar], T_{bp} [K] the boiling point under standard conditions, and *b* the solvent specific constant.

2.3 Evaporation of low-boiling solvents

2.3.1 Solvent: Acetone (≥99.8%)

| System Configura | ation | | | |
|---------------------------|---|---|--------------|-------------|
| Multivapor | P-12 | | | |
| Cover | PETP (black) | | | |
| Cover sealing | EPDM (black) | | | |
| Vessel type | BUCHI SpeedExtractor | | | |
| Options | Vacuum controller V-8 | | | |
| | Vacuum pump V-700 | | | |
| Parameters | | | | |
| Solvent | Acetone Evaporated volume Category Boiling point | ≥99.8% 240 ml, 20 ml/position Low-boiling 56°C | | |
| | Desired boiling point | 30°C | | |
| Heating medium | Temperature Volume | 55°C 45 ml H₂O/position | | |
| Cooling medium (F-105) | Water/glycol Temperature | 70:30% (v/v) 10°C | | |
| (1 100) | Measured flow | 1200 ml/min | Evaporation | performance |
| Rotational speed | Position 9 | | Positions | . 12 |
| | | | Total volume | 240 ml |
| Pressure gradien | t | | Time | 11 min |
| Starting point | 850 mbar | | Rate overall | 1.3 l/h |
| Ramp | 850-518 mbar | 1 min | Rate per | |
| | 518-361 mbar | 2 min | position | 109 ml/h |
| | 361-333 mbar | 2 min | pooldon | |
| Constant | 333 mbar | 6 min | | |
| _ | | | | |
| | 1000 | | | |



<u>S min</u>

Results

| ncound | |
|---------------------------------------|--------------|
| Total process time | 11 min |
| Solvent recovery after main condenser | 235 ml (98%) |
| Solvent recovery after secondary | < 1 ml |
| condenser | |
| | |

| EasyVac | Detection start point | Distillation | Detection end point |
|-----------------------------------|-----------------------------|--------------|----------------------------------|
| Acetone | \square | \square | $\mathbf{\overline{\mathbf{A}}}$ |
| Limitations: | | | |
| Glassware diameter < 13 mm | | | |
| Glassware diameter 13 < x < 19 mm | \square | \square | |
| Glassware diameter > 19 mm | \square | M | V |

2.3.2 Solvent: Ethyl acetate (Analytical grade)

| | System Configura | ation | | | |
|--|------------------|--|---|---------------|------------|
| Cover sealing Vessel type FFKM (black) BUCHI SpeedExtractor tube (4953) Options Vacuum controller V-855 (firmware 3.03) Vacuum pump V-700 Parameters Solvent Ethyl acetate Evaporated volume Category Desired boiling point Category Desired boiling point Category Desired boiling point Category Desired boiling point Colong medium (F-105) Ethyl acetate Evaporated volume S5°C Volume Analytical grade S5°C Volume Cooling medium (F-105) Temperature Measured flow 55°C 1200 ml/min Position Pressure gradient Total volume 242 mbar 2 min 2 min 242 mbar Evaporation performan Positions Evaporation performan Positions Constant 140 mbar 5 min Simin Results Total process time Solvent recovery after secondary 3 min 2 min 2 min Solvent recovery after secondary 10 min 2 min 2 min | Multivapor | P-12 | | | |
| Vessel type BUCHI SpeedExtractor tube (49535) Options Vacuum controller V-855 (firmware 3.03) Vacuum pump V-700 Parameters Solvent Ethyl acetate Evaporated volume 240 ml, 20 ml/position Category Boiling point Tomperature Analytical grade 240 ml, 20 ml/position Category Boiling point Tomperature First Station Cooling medium Vater/glycol First Station Cooling medium Water/glycol First Station Cooling medium Vater/glycol Evaporation before the assured flow Evaporation performan Rotational speed Position 9 Total volume 2 min 242 - 140 mbar 2 min 242 - 140 mbar Time 10 min 200 Constant 140 mbar 5 min 5 min 120 ml Results Total process time Solvent recovery after secondary 2 min 2 37 ml (98.8%) Solvent recovery after secondary 10 min 2 37 ml (98.8%) | | | | | |
| Options Vacuum controller V-855 (firmware 3.03) Vacuum pump V-700 Parameters Solvent Ethyl acetate Evaporated volume Category Analytical grade 240 ml, 20 ml/position Desired boiling point 30°C Heating medium (F-105) Temperature Temperature 55°C Volume Total volume 45 ml H ₂ O/position Rotational speed Position 9 Position 9 Pressure gradient Total volume 242 mbar 2 min 242 mbar Constant 140 mbar 5 min Constant 140 mbar 5 min Results Total process time Solvent recovery after main condenser Total process time Solvent recovery after secondary 10 min 237 ml (98.8%) Solvent recovery after secondary | | | | | |
| Vacuum pump V-700 Parameters Solvent Ethyl acetate Evaporated volume Category Boiling point Analytical grade 240 ml, 20 ml/position Low-boiling Feating medium (F-105) Temperature Volume Vater/glycol 25°C 70:30% (v/v) Temperature 10°C Measured flow 200 ml/min Rotational speed Position 9 Position 9 Pressure gradient Starting point 242:140 mbar 2 min 242:140 mbar 1 min 242:140 mbar Constant 140 mbar 5 min Total volume 200 ml/min Starting point 242 mbar 5 min Results Results Total process time Solvent recovery after main condenser 237 ml (98.8%) 237 ml (98.8%) | | | | | |
| Parameters Solvent Ethyl acetate Evaporated volume Category Analytical grade 240 ml, 20 ml/position Desired boiling point Trice Trice Total volume Analytical grade 240 ml, 20 ml/position Heating medium Temperature 55°C Volume 45 ml H ₂ O/position Evaporation performan Cooling medium Water/glycol 70:30% (i/v) Total volume 240 ml Rotational speed Position 9 Position 9 Positions 12 Total volume 242 mbar 2 min 242 mbar 1 min 242 mbar 2 min 242 mbar 2 min 242 mbar 2 min 242 mbar 2 min 242 mbar 2 min 242 mbar 2 min 242 mbar 2 min Constant 140 mbar 5 min 3 min Results Total process time 10 min Solvent recovery after main condenser 237 ml (98.8%) Solvent recovery after secondary | Options | | 55 (firmware 3.03) | | |
| Solvent Ethyl acetate Evaporated volume Category Desired boling point Temperature (F-105) Analytical grade 240 ml, 20 ml/position 30°C Heating medium Cooling medium (F-105) Temperature Temperature 10°C Measured flow 1200 ml/min 30°C 70:30% (v/v) Temperature 10°C Measured flow 1200 ml/min Evaporation performant Positions Pressure gradient Total volume 242 mbar 2 min 242 mbar 2 min 242 mbar Constant 140 mbar 5 min Constant 140 mbar 5 min | | Vacuum pump V-700 | | _ | |
| Evaporated volume Category Boiling point Category Desired boiling point Temperature (F-105) 240 ml/ 20 ml/position To Coling medium (F-105) Temperature Temperature Temperature 10°C Measured flow Temperature 10°C Measured flow 1200 ml/min Fereperature Position 9 Evaporation performant Positions Pressure gradient Starting point Starting point 242 mbar 653 mbar 653 -242 mbar 2 min 2 442 mbar Total volume 2 42 mbar 240 ml 1 min 2 422-140 mbar Constant 140 mbar 5 min Rate overall 1.4 l/h Rate per position 1000 000000000000000000000000000000000 | Parameters | | | | |
| Heating medium (Cooling medium (F-105) Temperature Volume 55°C 45 ml H ₂ O/position Rotational speed Position 9 Position 10°C Measured flow 1200 ml/min Pressure gradient 553 mbar Rate overall 1.4 l/h Rate overall Starting point 653 mbar 2 min Rate overall 1.4 l/h Rate per Constant 140 mbar 5 min Constant 140 mbar 5 min Results Total process time 10 min Solvent recovery after main condenser 237 ml (98.8%) Solvent recovery after secondary < 1 ml | Solvent | Evaporated volume Category Boiling point | 240 ml, 20 ml/position Low-boiling 77°C | | |
| (F-105) Temperature Measured flow 10°C 1200 ml/min Rotational speed Position 9 Pressure gradient Position 2 Starting point 653 nbar 653-242 mbar 2 min 242 mbar 1 min 242 nbar 2 min Constant 140 mbar 5 min 5 min Results Total process time train condenser 10 min 237 ml (98.8%) Solvent recovery after secondary 10 min 237 ml (98.8%) | - | Temperature Volume | 45 ml H ₂ O/position | | |
| Measured flow 1200 ml/min Rotational speed Position 9 Pressure gradient Position 1 Starting point 653 mbar 2 min 242 mbar 1 min 242-140 mbar 2 min Constant 140 mbar 5 min 12 Total volume 240 ml Image: Constant 140 mbar 5 min 700 120 ml Results Total process time 10 min 201 ml 203 ml 203 ml Results Total process time 10 min 203 ml 98.8%) 203 ml 98.8% | | | | | |
| Rotational speed Position 9 Positions 12 Pressure gradient Total volume 240 ml Starting point 653 mbar 2 min Armon 242 mbar 1 min 242 mbar 1 min 242 mbar 242 mbar 2 min Constant 140 mbar 5 min Image: Constant 140 mbar 6 8 10 12 Image: Constant 10 min <td>(F-105)</td> <td></td> <td></td> <td>Evaporation p</td> <td>orformance</td> | (F-105) | | | Evaporation p | orformance |
| Total volume 240 ml Starting point 653 mbar 2 min Ramp 653-242 mbar 2 min Rate overall 1.4 l/h Ramp 242 mbar 1 min 242-140 mbar 2 min Rate overall 1.4 l/h Constant 140 mbar 2 min 5 min 120 ml 10 min Total volume 240 ml 0 242-140 mbar 2 min 120 ml 120 ml 0 0 2 4 6 8 10 12 0 0 2 4 6 8 10 12 14 Results Total process time 10 min 237 ml (98.8%) 237 ml (98.8%) 50/vent recovery after secondary 2 min | Rotational speed | | 1200 11/11/1 | | |
| Starting point 653 mbar Rate overall 1.4 l/h Ramp 242 mbar 1 min 242 mbar 2 min Constant 140 mbar 2 min 200 120 ml Operation time (min) Results Total process time 10 min Solvent recovery after main condenser 237 ml (98.8%) Solvent recovery after secondary < 1 ml | | | | | 240 ml |
| Starting point 653 mbar Rate overall 1.4 l/h Ramp 653-242 mbar 2 min Rate overall 1.4 l/h Ramp 242 mbar 1 min 242-140 mbar 2 min Constant 140 mbar 5 min 120 ml Image: Constant 140 mbar 12 ml 120 ml Image: Constant 140 mbar 140 ml 140 ml Image: Constant 140 mbar 12 ml 120 ml Image: Constant 140 mbar 10 ml 120 ml Image: Constant 10 min 10 ml 120 ml | Pressure gradien | t | | Time | 10 min |
| Ramp 653-242 mbar 2 min 242 mbar 1 min 242-140 mbar 2 min Rate per position Constant 140 mbar 5 min Image: Constant 140 mbar 5 min 5 min Results Total process time 10 min Solvent recovery after main condenser 237 ml (98.8%) Solvent recovery after secondary < 1 ml | - | | | Rate overall | 1.4 l/h |
| 242-140 mbar 2 min Constant 140 mbar 5 min Image: Constant of the second and the second a | | 653-242 mbar | | Rate per | |
| Results Total process time 10 min Solvent recovery after main condenser 237 ml (98.8%) Solvent recovery after secondary < 1 ml | | | | position | 120 111/11 |
| Results Total process time 10 min Solvent recovery after main condenser 237 ml (98.8%) Solvent recovery after secondary < 1 ml | Ormationt | 4.40 mb an | E min | | |
| Results Total process time 10 min Solvent recovery after main condenser 237 ml (98.8%) Solvent recovery after secondary < 1 ml | Constant | 140 mbar | 5 min | | |
| Total process time10 minSolvent recovery after main condenser237 ml (98.8%)Solvent recovery after secondary< 1 ml | Lashin Disson | 600 500 400 200 100 0 2 4 6 | | | |
| Solvent recovery after main condenser237 ml (98.8%)Solvent recovery after secondary< 1 ml | | | | | |
| Solvent recovery after secondary < 1 ml | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

| EasyVac | Detection start | | Detection |
|-----------------------------------|----------------------------------|--------------|-----------|
| | point | Distillation | end point |
| Ethyl acetate | $\mathbf{\overline{\mathbf{A}}}$ | \square | \square |
| Limitations: | | | |
| Glassware diameter < 13 mm | | | |
| Glassware diameter 13 < x < 19 mm | \checkmark | \checkmark | |
| Glassware diameter > 19 mm | $\mathbf{\overline{\mathbf{A}}}$ | \square | \square |

2.3.3 Solvent: *n*-Hexane (purum ≥95.0%)

| System Configura | ation | | | |
|------------------|------------------------|---|---------------------------------------|--------------|
| Multivapor | P-12 | | | |
| Cover | PETP (black) | | | |
| Cover sealing | FKM (green) | | | |
| Vessel type | BUCHI SpeedExtractor | | | |
| Options | Vacuum controller V-85 | 55 (firmware 3.03) | | |
| | Vacuum pump V-700 | | _ | |
| Parameters | | | | |
| Solvent | n-Hexane | Purum ≥95.0% | | |
| | Evaporated volume | 360 ml, 30 ml/position | | |
| | Category | Low-boiling | | |
| | Boiling point | 69°C | | |
| | Desired boiling point | 30°C | | |
| Heating medium | Temperature | 55°C | | |
| Cooling medium | Volume Water/glycol | 45 ml H ₂ O/position 70:30% (v/v) | | |
| (F-105) | Temperature | 10°C | | |
| (1-103) | Measured flow | 1200 ml/min | Evaporation pe | erformance |
| Rotational speed | Position 9 | 1200 110/1111 | Positions | 12 |
| | | | Total volume | 1∠ 360 ml |
| Pressure gradien | t | | Time | 10 min |
| Starting point | 741 mbar | | Rate overall | 2.2 l/h |
| Ramp | 741-291 mbar | 2 min | Rate per | |
| | 291-237 mbar | 2 min | · · · · · · · · · · · · · · · · · · · | 180 ml/h |
| | 237-230 mbar | 2 min | position | |
| | 230-202 mbar | 3 min | | |
| Constant | 202 mbar | 1 min | | |
| Pressure [mharl | 0 | | | |

| Results | |
|---------------------------------------|---------------|
| Total process time | 10 min |
| Solvent recovery after main condenser | 360 ml (100%) |
| Solvent recovery after secondary | < 1 ml |
| condenser | |

| EasyVac | Detection start point | Distillation | Detection end point |
|-----------------------------------|-----------------------------|-------------------------|------------------------|
| <i>n</i> -Hexane | | $\mathbf{\overline{A}}$ | \checkmark |
| Limitations: | | | |
| Glassware diameter < 13 mm | | | |
| Glassware diameter 13 < x < 19 mm | \checkmark | \checkmark | |
| Glassware diameter > 19 mm | \square | \square | |

2.3.4 Solvent: *n*-Hexane – P-6 (purum ≥95.0%)

| System Configu | ration | | | |
|---------------------------|-----------------------------|-------------------------|---------------|------------|
| Multivapor | P-6 | | | |
| Cover | PETP (black) | | | |
| Cover sealing | FKM (green) | | | |
| Vessel type | BUCHI SpeedExtractor | | | |
| Options | Vacuum controller V-85 | 55 (firmware 3.03) | | |
| | Vacuum pump V-700 | | | |
| | | | | |
| Parameters | | | | |
| Solvent | <i>n</i> -Hexane | Purum ≥95.0% | | |
| | Evaporated volume | 600 ml, 100 ml/position | | |
| | Category | Low-boiling | | |
| | Boiling point | 69°C | | |
| | Desired boiling point | 30°C | | |
| Heating medium | Temperature | 55°C | | |
| | Volume | 145 ml H_2O /position | | |
| Cooling medium (F-105) | Water/glycol Temperature | 70:30% (v/v) 10°C | | |
| (F-105) | Measured flow | 1200 ml/min | Evaporation p | erformance |
| Rotational speed | | | Positions | 6 |
| | | | Total volume | 600 ml |
| Pressure gradie | nt | | Time | 14 min |
| _ | 741 mbar | | Rate overall | 2.6 l/h |
| Starting point Ramp | 741 mbar 741-308 mbar | 2 min | | 2.0 1/11 |
| ιταπρ | 308-241 mbar | 2 min | Rate per | 429 ml/h |
| | 241 mbar | 3 min | position | 425 111/11 |
| | 241-213 mbar | 1 min | | |
| Constant | 213 mbar | 6 min | | |
| | | | | |
| | | | | |
| | 800 - | | | |
| | F ⁷⁰⁰ | | | |
| | Lengua 200 | | | |
| | e 400 | | | |
| | 300 | | | |
| | e 200 | | | |
| | | | | |
| | | 10 15 | | |
| | Evaporation time | [min] | | |
| L | | | | |
| | | | | |
| | | | | |

| Total process time | 14 min |
|---------------------------------------|----------------|
| Solvent recovery after main condenser | 583 ml (97.2%) |
| Solvent recovery after secondary | 2 ml |
| condenser | |

| EasyVac | Detection | | |
|------------------|-----------|--------------|-----------|
| | start | | Detection |
| | point | Distillation | end point |
| <i>n</i> -Hexane | \square | (☑) | Too early |

Evaporation of medium-boiling solvents 2.4

| 2.4.1 | Solvent: X | ylene (| Mixture of isc | omers, purum ≥97.0%) | |
|-------|------------|---------|----------------|----------------------|--|
|-------|------------|---------|----------------|----------------------|--|

| | ition | | | | | |
|---|--|---|-----------------|------------------------|-------------------------------|------------|
| Multivapor | P-12 | | | | | |
| Cover | PETP (black) | | | | | |
| Cover sealing | FKM (green) | | | | | |
| Vessel type | BUCHI SpeedExtrac | tor tube (49535 | 5) | | | |
| Options | Vacuum controller V- | | | | | |
| | Vacuum pump V-700 |) | | | | |
| | Woulff flask (48857) | | | | | |
| Parameters | | | | | | |
| Solvent | Xylene | Purum ≥9 | 7 0º/ (icor | more) | | |
| Solveni | | 240 ml, 20 | | | | |
| | Evaporated volume | Medium-b | | UII | | |
| | Category | 139°C | oning | | | |
| | Boiling point | | | | | |
| L la attinar un a divun | Desired boiling point | 50°C | | | | |
| Heating medium | Temperature | 80°C | \/= = = :+! = : | | | |
| Caaling and l | Volume | 45 ml H ₂ C | | | | |
| Cooling medium | Water/glycol | 70:30% (v | //V) | | | |
| (F-105) | Temperature | 10°C | | | F orman and the second | |
| | Measured flow | 1200 ml/m | าเท | | Evaporation p | erformance |
| Rotational speed | Position 9 | | | | Positions | 12 |
| | | | | | Total volume | 240 ml |
| Pressure gradient | t | | | | Time | 7 min |
| Starting point | 200 mbar | | | | Rate overall | 2.1 l/h |
| Ramp | 200-30 mbar | 3 min | | | Rate per | 474 10 |
| | | | | | position | 171 ml/h |
| Constant | 30 mbar | 4 min | | | | |
| Constant | 30 mbar 250 150 100 50 0 2 4 Evaporation ti | 4 min | | | | |
| Results | 250 200 150 50 0 2 4 Evaporation ti | 6 8 me [min] | | | | |
| Results Total process time | 250 200 150 50 0 2 4 Evaporation ti | 6 8 me [min] 7 min | | | | |
| Results Total process time Solvent recovery a | 250 150 100 50 0 2 4 Evaporation ti fter main condenser | 6 8 me [min] 7 min 187 ml (77.9% |) | | | |
| Results Total process time Solvent recovery a Solvent recovery in | 250 150 100 50 0 2 4 Evaporation ti fter main condenser Woulff flask | 6 8 me [min] 7 min 187 ml (77.9% 52 ml (21.7%) |) | | | |
| Results Total process time Solvent recovery a Solvent recovery a Solvent recovery a | 250 150 100 50 0 2 4 Evaporation ti fter main condenser Woulff flask | 6 8 me [min] 7 min 187 ml (77.9% |) | | | |
| Results Total process time Solvent recovery a Solvent recovery in | 250 150 100 50 0 2 4 Evaporation ti fter main condenser Woulff flask | 6 8 me [min] 7 min 187 ml (77.9% 52 ml (21.7%) |) | | | |
| Results Total process time Solvent recovery an Solvent recovery an Solvent recovery an condenser | 250 150 100 50 0 2 4 Evaporation ti fter main condenser Woulff flask | 6 8 me [min] 7 min 187 ml (77.9% 52 ml (21.7%) |) | | | |
| Results Total process time Solvent recovery a Solvent recovery a Solvent recovery a | 250 150 100 50 0 2 4 Evaporation ti fter main condenser Woulff flask | 7 min 187 ml (77.9% 52 ml (21.7%) 1 ml |) | Detection | | |
| Results Total process time Solvent recovery an Solvent recovery an Solvent recovery an condenser | 250 150 100 50 0 2 4 Evaporation ti fter main condenser Woulff flask | 7 min 187 ml (77.9% 52 ml (21.7%) 1 ml Detection start |) stillation | Detection end point | | |

Evaporation of high-boiling solvents 2.5

| Multivapor | P-12 | | | |
|------------------|---|---|----------------|---------|
| Cover | PETP (black) | | | |
| Cover sealing | EPDM (black) | | | |
| Vessel type | BUCHI SpeedExtractor | | | |
| Options | Vacuum controller V-85 | 55 (firmware 3.03) | | |
| | Vacuum pump V-700 | | | |
| Parameters | | | | |
| Solvent | N,N-Dimethyl- formamide (DMF) Evaporated volume Category Boiling point Desired boiling point | Purum ≥98.0% 240 ml, 20 ml/position High-boiling 153°C 60°C | | |
| Heating medium | Temperature Volume | 80°C 45 ml H₂O/position | | |
| Cooling medium | Water/glycol | 70:30% (v/v) | | |
| (F-105) | Temperature | 20°C | Evaporation pe | erforma |
| | Measured flow | 1200 ml/min | Positions | 12 |
| Rotational speed | Position 10 | | Total volume | 240 |
| | | | Time | 25 m |
| Pressure gradien | t | | Rate overall | 0.6 l |
| Starting point | 400 mbar | | Rate per | |
| Ramp | 400-28 mbar | 2 min | position | 48 m |
| | | | | |
| Constant | 28 mbar | 23 min | _ | |
| Γ | 500 - | | _ | |
| | 400 300 200 100 0 | | | |
| | | 20 30 | | |
| | · · · · | | | |

Solvent: Dimethylformamide, DMF (purum, ≥98.0%) 2.5.1

| Results |
|---------|
|---------|


| Total process time | 25 min |
|---------------------------------------|----------------|
| Solvent recovery after main condenser | 236 ml (98.3%) |
| Solvent recovery after secondary | 1 ml |
| condenser | |

| EasyVac | Detection start point | Distillation | Detection end point |
|-----------------------|-----------------------------|--------------|------------------------|
| N,N-Dimethylformamide | | | |

2.6 Evaporation of polar protic solvents

2.6.1 Solvent: Ethanol (4.8% methanol, purum, ≥99.8%)

| System Configura | ation | |
|-------------------|--------------------------|---------------------------------|
| Multivapor | P-12 | |
| Cover | PETP (black) | |
| Cover sealing | EPDM (black) / FKM (gree | en) |
| Vessel type | BUCHI SpeedExtractor tu | |
| Options | Vacuum controller V-855 | (firmware 3.03) |
| | Vacuum pump V-700 | |
| | PE frit (44856) | |
| | | |
| Parameters | | |
| Solvent | Ethanol (4.8% methanol) | Purum, ≥99.8% |
| | Evaporated volume | 240 ml, 20 ml/position |
| | Category | Polar protic |
| | Boiling point | 78°C |
| | Desired boiling point | 30°C |
| Heating medium | Temperature | 55°C |
| | Volume | 45 ml H ₂ O/position |
| Cooling medium | Water/glycol | 70:30% (v/v) |
| (F-105) | Temperature | 10°C |
| | Measured flow | 1200 ml/min |
| Rotational speed | Position 9 | |
| | | |
| Pressure gradient | t | |
| Starting point | 602 mbar | |
| Ramp | 602-173 mbar | 2 min |
| | 173-96 mbar | 3 min |
| | 96-86 mbar | 4 min |
| | 86 mbar | 4 min |
| | 86-78 mbar | 6 min |



| 12 |
|---------|
| 240 ml |
| 19 min |
| 0.8 l/h |
| 63 ml/h |
| |

Evaporation performance



| Total process time | 19 min |
|---------------------------------------|----------------|
| Solvent recovery after main condenser | 239 ml (99.6%) |
| Solvent recovery after secondary | < 1 ml |
| condenser | |

| EasyVac | Detection start point | Distillation | Detection end point |
|-----------------------------------|----------------------------------|--------------|------------------------|
| Ethanol | \square | \square | \square |
| Limitations: | | | |
| Glassware diameter 13 < x < 19 mm | $\mathbf{\overline{\mathbf{A}}}$ | \square | |
| Glassware diameter > 19 mm | \checkmark | \square | \square |

2.6.2 Solvent: Distilled water, H₂O

| System Configura | tion | | | | | |
|---------------------|--|------------------|---------------------------|-----------|----------------|------------|
| Multivapor | P-12 | | | | | |
| Cover | PETP (black) | | | | | |
| Cover sealing | EPDM (black) / FKN | Л (green) | | | | |
| Vessel type | BUCHI SpeedExtra | | | | | |
| Options | Vacuum controller \ | /-855 (firmwa | are 3.03) | | | |
| | Vacuum pump V-70 | | | | | |
| | Woulff flask (48857 |) | | | | |
| | | | | | | |
| Parameters | Distillations for | DUOU | | 04.0 | | |
| Solvent | Distilled water | | I Fontavapor | | | |
| | Evaporated volume | | I, 20 ml/positi | ion | | |
| | Category | Polar J | | | | |
| | Boiling point | 100°C | | | | |
| | Desired boiling poin | | | | | |
| Heating medium | Temperature | 55°C | | | | |
| | Volume | | H ₂ O/position | | | |
| Cooling medium | Water/glycol | 70:309 | % (v/v) | | | |
| (F-105) | Temperature | 10°C | | | | |
| | Measured flow | 1200 r | nl/min | | Evaporation pe | erformance |
| Rotational speed | Position 9 | | | | Positions | 12 |
| | | | | | Total volume | 240 ml |
| Pressure gradient | | | | | Time | 67 min |
| Starting point | 542 mbar | | | | Rate overall | 0.2 l/h |
| Ramp | 542-69 mbar | 2 mir | | | Rate per | |
| _ | 69-33 mbar | 3 mir | | | position | 18 ml/h |
| Constant | 33 mbar | 57 m | | | position | |
| Constant | 33-27 mbar 27 mbar | 2 mir 3 mir | | | | |
| Pressure [mbar] | 600 500 400 300 200 100 0 20 40 Evaporation | 60 time [min] | 80 | | | |
| Deculto | | | | | | |
| Results | | 07 . | | _ | | |
| Total process time | tor main sand- | 67 min | 00() | | | |
| | ter main condenser | 235 ml (97 | | | | |
| Solvent recovery in | | 3 ml (1.3%) |) | | | |
| Solvent recovery at | ter secondary | 2 ml | | | | |
| condenser | | | | | | |
| | | | | | | |
| EasyVac | | Detection | | Dotootion | | |
| | | start | Distill | Detection | | |
| | | point | Distillation | end point | | |
| Water | | V | | | | |
| Limitations: | | | | | | |
| Glassware diamete | r < 13 mm | | | | | |
| | - 10 · · · · · 10 · | | | | | |
| Glassware diamete | | র্ ব | | | | |

3 Evaporation of solvent mixtures

Evaporation of the solvent mixtures hexane / dichlormethane and ethanol / water has been discussed extensively in the best@buchi 06/2001.⁷ Custom mixtures are addressed similarly. Solvent mixtures are usually distilled by evaporation of one component after the other. In the presented examples, a boiling temperature of 30 to 40°C is sought in order to protect the products.

Optimal conditions are determined in a boiling diagram. Theoretical useful pressures are determined by the heat medium temperature and the lowest possible condensation temperature. The boiling curve is calculated using Equation 1 (Section 2.2).

3.1 Hexane / dichloromethane 7:3 (e.g. flash chromatography)

The problem lies in the fact that the dichloromethane has to be distilled off at a relatively high system pressure of from 750 to 800 mbar, while hexane requires a system pressure of about 300 mbar (**Figure 10**). If dichloromethane is distilled off at too low pressure, then the condenser will soon be overburdened and solvent emissions into the environment become unavoidable.

To avoid this, dichloromethane is distilled off at an only slightly reduced pressure and then the system pressure is lowered to distill hexane. Such a lowering of pressure in stages can of course be done manually.

The monitoring time is considerable, however. Due to the fact that the boiling temperature of the hexane residue rises at a certain pressure with decreasing dichloromethane content, there is also always the danger that there could be a boiling delay if a stage is slightly too large.

If, however, the system pressure is continuously lowered after a short starting period, i.e., along a pressure gradient, then this danger is avoided and the distillation process can be carried out without supervision or servicing.




Figure 10: Boiling diagram for dichloromethane/hexane. *a* Heat medium temperature; *b* Lowest possible condensation temperature for a cooling temperature of 10°C. The rectangle area, which is surrounded in green, marks the theoretically useful area; $c \rightarrow d$ shows the used pressure range.

⁷ A. Talamona, G. Schönenberger, Vacuum control during vaporization of problematic solutions, best@buchi no. 06/2001.

Solvent mixture: Hexane / dichloromethane 7:3

| | ation | | | |
|----------------------------|----------------------------------|---------------------------------|-------------------|-----------|
| Multivapor | P-12 | | | |
| Cover | PEEK (ecru) / PETP (bl | lack) | | |
| Cover sealing | FKM (green) | | | |
| Vessel type | BUCHI SpeedExtractor | | | |
| Options | Vacuum controller V-85 | 55 (firmware 3.03) | | |
| | Vacuum pump V-700 | | | |
| Parameters | | | | |
| Solvent mixture | Hexane | ≥96% | | |
| | Dichloromethane | ≥99.8% | | |
| | Ratio | 7:3 (v/v) | | |
| Heating medium | Evaporated volume Temperature | 240 ml, 20 ml/position 55°C | | |
| nealing meulum | Volume | 45 ml H ₂ O/position | | |
| Cooling medium | Water/glycol | 70:30% (v/v) | | |
| (F-105) | Temperature | 10°C | | |
| (1 100) | Measured flow | 1200 ml/min | | |
| Rotational speed | Position 9 | | Evaporation pe | erformanc |
| | | | Positions | 12 |
| Pressure gradien | t | | Total volume | 240 ml |
| 0, | 760 mbar | | Time | 12 min |
| Starting point | 760 mbar | 3 min | Rate overall | 1.2 l/h |
| Starting point Constant | 700 mbai | | | |
| | 760-350 mbar | 7 min | Rate per | 100 ml/ |
| Constant | | | Rate per position | 100 ml/ |

| 0 | 3 | 6 | 9 | 12 | 15 |
|---|-----|----------|-----------|------|----|
| | Eva | aporatio | n time [ı | nin] | |

| Results | |
|---------------------------------------|----------------|
| Total process time | 12 min |
| Solvent recovery after main condenser | 237 ml (98.8%) |
| Solvent recovery after secondary | 1 ml |
| condenser | |

3.2 Ethanol / water 6:4 (e.g. reverse phase-chromatography)

Ethanol/water mixtures tend to bump and foam during the distillation process. If the distillation process is carried out very carefully, then the process can be almost endless. If, on the other hand, the pressure is lowered a bit too sharply, then the solution foams over and the distillation has to be started all over again! Additionally, the distillation is performed at a truly low-pressure level to achieve the desired boiling point (**Figure 11**).

In this application as well, a vacuum control system with programmable pressure gradient is very helpful and results in a significant reduction in laboratory personnel. The distillation record sheet on the next page shows how the pressure during manual operation is reduced rather tediously during the first 20 minutes of operation, always close to the limit of over-foaming, and how the problem can be elegantly solved with the use of a gradient function. True, the distillation process lasts a bit longer, but the danger of over-foaming is clearly eliminated.





Figure 11: Boiling diagram for ethanol/water. *a* Heat medium temperature; *b* Lowest possible condensation temperature for a cooling temperature of 10°C. The rectangle area, which is surrounded in green, marks the theoretically useful area; $c \rightarrow d$ shows the used pressure range.

Solvent mixture: Ethanol / water 6:4

| System Configur | ation | |
|-----------------|---|-----------------------|
| Multivapor | P-12 | |
| Cover | PETP (black) | |
| Cover sealing | EPDM (black) / FKM (gre | en) |
| Vessel type | BUCHI SpeedExtractor tube (49535) | |
| Options | Vacuum controller V-855 (firmware 3.03) | |
| | Vacuum pump V-700 | |
| | | |
| Parameters | | |
| Solvent mixture | Ethanol (absolute) | purum ≥99.8% |
| | Distilled water | BUCHI Fontavapor 210 |
| | Ratio | 6:4 (v/v) |
| | Evanorated volume | 240 ml 20 ml/position |

| Biotinou nator | Been in ontarapor 210 |
|-------------------|---|
| Ratio | 6:4 (v/v) |
| Evaporated volume | 240 ml, 20 ml/position |
| Temperature | 55°C |
| Volume | 45 ml H₂O/position |
| Water/glycol | 70:30% (v/v) |
| Temperature | 10°C |
| Measured flow | 1200 ml/min |
| Position 9 | |
| | Ratio Evaporated volume Temperature Volume Water/glycol Temperature Measured flow |

| Pressure gradient | | | |
|-------------------|--------------|--------|--|
| Starting point | 602 mbar | | |
| Ramp | 602-267 mbar | 2 min | |
| | 267-107 mbar | 4 min | |
| | 107-84 mbar | 2 min | |
| | | | |
| Constant | 84 mbar | 49 min | |



| Evaporation performance | | |
|-------------------------|---------|--|
| Positions | 12 | |
| Total volume | 240 ml | |
| Time | 57 min | |
| Rate overall | 0.3 l/h | |
| Rate per position | 21 ml/h | |

| Results | |
|---------|--|
| | |

| Total process time | 57 min |
|---------------------------------------|----------------|
| Solvent recovery after main condenser | 225 ml (93.8%) |
| Solvent recovery after secondary | 4 ml |
| condenser | |

3.3 Azeotrope chloroform / methanol 8:1 (e.g. flash chromatography)

An azeotrope a mixture of two or more liquids in such a ratio that its composition cannot be changed by simple distillation. The vapour of boiled azeotrope has the same ratio of constituents as the original mixture

Each azeotrope has a characteristic boiling point. The boiling point temperature of an azeotrope is either less than the boiling point temperatures of any of its constituents (a positive azeotrope), or greater than the boiling point temperatures of any of its constituents (a negative azeotrope). Decrease of boiling point is seen in the phase diagrams of the positive azeotrope chlororform / methanol (**Figure 12**).




Figure 12: Phase diagrams of the positive azeotrope chlororform / methanol.⁸

If two solvents can form a positive azeotrope, then distillation of any mixture of those constituents will result in the distillate being closer in composition to the azeotrope than the starting mixture. Repeated series distillations will produce mixtures that are progressively closer to the azeotropic ratio. But you cannot concentrate above the azeotropic ratio. Distillation of a mixture above the azeotropic point, the distillate (contrary to intuition) will be poorer in higher concentrated volatile compound than the original but slightly richer than the azeotrope.

Similarly, if two solvents can form a negative azeotrope, then distillation of any mixture of those constituents will result in the residue being closer in composition to the azeotrope than the original mixture. Boiling of any negative azeotrope solution long enough will cause the solution left behind to approach the azeotropic ratio.

Non-azeotropic solvent mixtures are distilled by evaporation of one component after the other. Assuming chloroform / methanol is non-azeotropic, optimal conditions are displayed in the rectangle area of the boiling diagram (**Figure 13**). In fact, chloroform / methanol at a ratio of

⁸ Azeotrope. Wikipedia.org, the free encyclopedia. http://en.wikipedia.org/wiki/Azeotrope. Last visit: February 17th 2011.

8:1 corresponding to 0.8 mol-% has a boiling point of 54°C (**Figure 12**). Therefore, distillation is performed at a higher vacuum than with non-azeotropic solvent mixtures ($c \rightarrow d$, **Figure 13**).





Figure 13: Boiling diagram for chloroform/methanol. *a* Heat medium temperature; *b* Lowest possible condensation temperature for a cooling temperature of 5°C. The rectangle area, which is surrounded in green, marks the theoretically useful area; $c \rightarrow d$ shows the used pressure range.

Solvent mixture: Chloroform / methanol 8:1

| System Configura | tion | | _ | |
|-------------------|-----------------------------|------------------------|---------------|--------|
| Multivapor | P-12 | | | |
| Cover | PEEK (ecru) / PETP (black) | | | |
| Cover sealing | FKM (green) / FFKM (bla | ack) | | |
| Vessel type | BUCHI SpeedExtractor | | | |
| Options | Vacuum controller V-858 | 5 (firmware 3.03) | | |
| | Vacuum pump V-700 | | | |
| | | | | |
| Parameters | | | _ | |
| Solvent mixture | Chloroform | puriss. p.a. ≥99.8% | | |
| | Methanol | ≥99.8% | | |
| | Ratio | 8:1 (v/v) | | |
| | Evaporated volume | 240 ml, 20 ml/position | | |
| Heating medium | Temperature | 55°C | | |
| | Volume | 45 ml H_2O /position | | |
| Cooling medium | Water/glycol | 70:30% (v/v) | | |
| (F-105) | Temperature | 5°C | | |
| Detetional anead | Measured flow Position 9 | 1200 ml/min | Evaporation p | orform |
| Rotational speed | Position 9 | | | |
| | | | Positions | 12 |
| Pressure gradient | t | | Total volume | 240 |
| Starting point | 800 mbar | | Time | 14 m |
| Ramp | 800-612 mbar | 1 min | Rate overall | 1.0 l |
| | 612-480 mbar | 2 min | Rate per | |
| | 480-400 mbar | 5 min | position | 86 m |
| Constant | 400 mbar | 6 min | | |
| Constant | | U IIIII | _ | |



| Evaporation performance | | |
|-------------------------|--|--|
| 12 | | |
| 240 ml | | |
| 14 min | | |
| 1.0 l/h | | |
| 86 ml/h | | |
| | | |

Results

| Total process time | 14 min |
|---------------------------------------|----------------|
| Solvent recovery after main condenser | 237 ml (98.8%) |
| Solvent recovery after secondary | 1 ml |
| condenser | |

4 Appendix A – Reference BUCHI Application Notes

4.1 Fat Determination in Shortbread using SpeedExtractor E-916

The determination of fat in food and feed is a routine procedure for quality assurance and labeling. Fat was determined in butter shortbread after homogenization of the sample and extraction with the SpeedExtractor E-916. The total fat content was determined gravimetrically after the extract has been dried to a constant weight. The fat content of 10.21% corresponds to the value determined with classic Soxhlet.


Introduction


The determination of fat in food and feed is a routine procedure for in quality assurance and labeling. Fat was determined in butter shortbread after homogenization of the sample and extraction with the SpeedExtractor E-916. The solvent was evaporated in parallel using the Multivapor P-6. The total fat content was determined gravimetrically after the extract has been dried to a constant weight.



Figure 14: Shortbread

| Extraction settings | | |
|---------------------|---|--|
| Temperature | 100°C | |
| Pressure | 100 bar | |
| Solvent | <i>n</i> -Hexane | |
| Cells | 40 ml | |
| Vials | 240 ml | |
| Cycles | 3 | |
| Heat-up | 1 min | |
| Hold | 5 min | |
| Discharge | 4/3/3 min (1 st /2 nd /3 rd cycle) | |
| Flush with solvent | 2 min | |
| Flush with gas | 4 min | |





Experimental

Instrumentation: Mixer B-400, SpeedExtractor E-916, Multivapor P-6 with Vacuum pump V-700 and Controller V-855, drying oven. The homogenized sample was mixed with quartz sand and extracted with the SpeedExtractor using the parameters shown in the adjacent table. The sample was extracted in triplicate.

The total time for the extraction is approx. 60 min and per position, approx. 50 ml solvent are used. The solvent was evaporated in parallel using the Multivapor P-6 (see **Figure 15**). The extracts were then dried to a constant weight in a drying oven (102°C) and the fat content was calculated.

Results

The determined fat content (**Table 1**) of 10.20% corresponds to the content obtained when extracting the same sample with Soxhlet extraction. With a fat content of 10.20% (rsd=0.42%, n=3) was determined.

 Table 1: Determined fat contents in short bread

| Results | Fat content |
|------------|-------------|
| Sample 1 | 10.25 % |
| Sample 2 | 10.14 % |
| Sample 3 | 10.24 % |
| Mean value | 10.21 % |
| rsd% | 0.60 |
| | |

References

SpeedExtractor E-916 operation manual For more detailed information refer to Application note 005/2009

4.2 Flash chromatography on cartridges for the separation of plant extracts

Flash chromatography on cartridges for the separation of plant extracts - Rules for the selection of chromatographic conditions, and comparison with MPLC

Weber P¹, Schafroth N², Hamburger M¹, Potterat O¹

¹ Division of Pharmaceutical Biology, University of Basel, Klingelbergstrasse, CH-4056 Basel, Switzerland

² Büchi Labortechnik AG, CH-9230 Flawil, Switzerland

Introduction

During the last decade, several systems for rapid preparative chromatography with pre-packed cartridges have been commercialized. Pre-packed cartridges ensure rapid separation cycles and ease of use.


separation cycles and ease of use. Flash chromotography systems were initially developed for rapid and easy purification of synthetic products, and numerous applications are documented in the experimental section of publications in synthetic chemistry. In contrast, application of such systems in the separation of complex natural product mixtures such as extracts has been nealected

The aim of this work was to explore the potential and limitations of cartridges for the In a wind that work was to explore the potential and minimums of candidges of the purification of natural product extracts. Empirical rules have been established for the determination of the separation conditions by preliminary TLC and HPLC analyses. The influence of solid introduction compared to liquid injection has been also investigated. The performance of the cartridges was compared to that of classical selfpacked MPLC glass columns for the separation of complex plant extracts of medicinal importance

Results

Correlation between HPLC and reversed phase flash chromatography

In preliminary experiments with reference compounds, a correlation was found between the capacity factors k' on the HPLC column and Sepacore cartridges $(K_{\text{Sepacore}} = 1.15 \times K_{\text{int},\text{C}}; R^2 = 0.97)$. HPLC separations can be transposed by increasing the gradient time by a factor 2-4. Theses rules have been applied to the separation of dicloromethane extracts of *Curcuma xanthorrhiza* and *Piper nigrum*.







Fig 1. Separation of a DCM extract of C. xanthorrhiza on RP-18.

A) HPLC column (4 x 150 mm) B) Sepacore® cartridge 40 x 150 mm C) glass column 26 x 240 mm $_{\rm 2}$ gives sortum to x 240 mm Mobile phase: H_O (A), MeOH (B) both with 0.1% HCOH; gradient: 50 – 100% B; flow rate: A) 1m/min; B/C) 30 mL/min; D) 10m/min; samples: A) 100 μ g in 10 μ L, B/C) 500 mg in 5 mL; detection at 220 nm.





6 MR 4.5-DRyd 7 M7 ARR Poertine

Fig. 2. Separation of a DCM extract of *P. nigrum* on RP-18.

RP-18. A) HPLC column (4 x 150 mm) B/C) Sepacore® cartridge (40 x 150 mm) D) Glass column 26 x 240 mm Mobile phase: H₂O (A), MeOH (B) both with 0.1% HCOOH. Gradient: 50 – 100% B; flow rate: A) 1mU/min, B/C(3) on L/ min, D) 10mU/min, Samples: A) 100 μg in 10 μL, B/C/D) 500 mg in 5 mL; detection at 270 nm.

Transposition from TLC to flash chromatography on silica gel

In preliminary experiments, a correlation was found between $1/R_i$ and k' values (k'= 1.1 x $1/R_i$, R' = 0.90). Mobile phase compositions for starting and end point of a gradient should be selected such as to obtain R_i values of 0.15-0.2 for the most Ipophilic component and the most hydrophilic constituent, respectively. This principle has been applied to the separation of a dichloromethane extract of Salvia miltiorrhiza.




Fig. 3. Separation of a DCM extract of S. miltiorrhiza on silica gel.

r-rg. -s. osparaturun in a Com extract 01 5. *Imitorimiza* on suică geli. Al Sepacorea cartidage (40 x 150 mm), B) Glass column 26 x 240 mm. Mobile phase: A/B) n-hexane (A), EtOAc (B); gradient: 2-50% B: slow rate: A) 30 mL/min, B) 10 mL/min; sample: A) 500 mg in 10 mL *n*-hexane-EtOAc (64). B) solid introduction using a precolumn; detection at 255 nm. C) TLC analysis of the extract and the fractions obtained from B with *n*-hexane-EtOAc MeOH (25:25:1); R1: tanshinone IIA, R2: cryptotanshinone.

Influence of the mode of sample introduction




Fig. 4. Influence of the sample introduction mode on the separation of a DCM extract of *S. milthiorrhiza*.

A) HPLC column (4 x 150 mm) B/C) Sepacore® cartridge (40 x 150 mm)

BC) Sepacoree carringge (40 x 150 mm) Mobile phase: H₀ (A), MeOH (B) both with 0.1% HCOOH; gradient A) 60 – 75% B in 30 min, BC/ 60% B for 10 min, 60 – 80% B in 110 min; flow rate: A) 1 mL/min, BC/ 30 mL/mir; detection at 250 nm; sample: A) 100 µg in 10 µL, B) liquid introduction, 500 mg in 8 mL MeOH, C) solid introduction, 500 mg.

Conclusions

· Reversed phase HPLC separations can be transposed by increasing the gradient time by a factor 2-4.

For normal phase separations, solvent compositions resulting in R_r values of 0.15-0.2 on TLC for the most lipophilic and the most hydrophilic constituents, respectively, should be selected as gradient endpoints.

 Sepacore® cartridges enabled a good separation of compounds with a broad range
of polarity, as typically found in plant extracts. The chromatographic resolution
remained, however, lowever than that achieved by MPLC on columns packed with
material of smaller particle size. For poorly soluble extracts, solid introduction gave
material extracts and the second s better results than liquid injection.

 Despite lower resolution as compared to MPLC, pre-packed cartridges are an attractive alternative for the purification of extracts and crude fractions due to their ease of use and speed of separation.

Instrumentation

raphy separations were performed on pre-packed silica (12 x 150 or 40 x 150 mm, *Büchi*) at a flow rate of ges), resp. Medium pressure liquid chromatography (MF

References

The original Poster is available at: http://www.buchi.com/Applications.20162.0.html

Publication:

Weber, P., Hamburger, M., Schafroth, N., Potterat, O. 2011. Flash chromatography on cartridges for the separation of plant extracts: Rules for the selection of chromatographic conditions and comparison with medium pressure liquid chromatography. Fitoterapia 82:155–161.

5 Appendix B – Recommended Multivapor version for a specific solvent

5.1 Abbreviations

| EPDM | Ethylenepropylenedimonomer | O-ring |
|--------|----------------------------|--------------------------|
| FKM | Fluoric caoutchouc | O-ring |
| FFKM | Perfluoro caoutchouc | O-ring |
| PEEK | Polyetheretherketone | Vacuum cover and adapter |
| PET(P) | Polyethyletherphtalate | Vacuum cover and adapter |
| PFA | Perfluoroalkoxy | Vacuum hose |
| PTFE | Polytetrafluoroethylene | Sealing discs |

5.2 Chemical resistance table

| | EPDM* ⁹ | FKM* ¹⁰ | FKKM* ¹¹ | PEEK* ⁹ | PET(P)* ¹³ | PFA** ¹² | PTFE** ¹² |
|----------------------------|--------------------|--------------------|---------------------|--------------------|-----------------------|---------------------|----------------------|
| Acetaldehyde | В | D | А | А | А | А | А |
| Acetic Acid | А | В | А | А | А | А | А |
| Acetic acid anhydride | В | D | А | А | А | А | А |
| Acetone | А | D | А | А | В | А | А |
| Benzene | D | А | А | А | А | А | А |
| Butanol | В | А | А | А | В | А | А |
| Chloroform | D | А | А | А | В | А | А |
| Diethyl ether | С | С | А | А | А | А | А |
| Dimethylformamide | А | - | А | А | В | А | А |
| Dimethylbenzene (Xylol) | D | A | A | A | A | A | А |
| Dioxane | В | - | А | А | А | А | А |
| Ethanol | А | А | А | А | А | А | А |
| Ethyl acetate | В | D | А | А | - | А | А |
| Hexane | С | А | А | А | А | А | А |
| Isobutanol | А | А | А | А | А | А | А |
| Isopropanol | А | А | А | А | А | А | А |
| Methanol | А | В | А | А | А | А | А |
| Methylene chloride | D | А | А | А | D | А | А |
| Nitrobenzene | С | В | А | В | D | А | А |
| Phenol | В | А | А | В | С | А | А |
| Propanol | А | А | А | А | А | А | А |
| Sulphuric acid, fuming | С | А | А | С | С | А | А |
| Carbon tetrachloride | D | А | А | А | А | А | А |

⁹ Operation Manual, Syncore® Accessories

¹⁰ Connectors, Verbindungstechnik AG, Chem. Beständigkeitsliste

¹¹ ExpotechUSA, www.expotechusa.com

¹² Semadeni, Chemical resistance table for polymers, www.semadeni.com

¹³ Prelon GmbH Dichtungstechnik, Chem. Beständigkeit von Kunststoffen, www.prelon.de

| | EPDM | FKM | FKKM | PEEK | PET(P) | PFA | PTFE |
|----------------------|------|-------|------|------|-----------------|------|------|
| | | I INN | | | F L 1(F) | FI A | FILE |
| Tetrahydrofurane | В | D | А | А | А | А | А |
| Toluene | D | А | А | А | А | А | А |
| Triethylamine | С | А | А | А | - | - | - |
| Trichloroethane | D | А | А | А | А | - | - |
| Trichloroacetic acid | В | - | - | А | - | А | А |
| Vinylidene chloride | D | - | А | А | В | - | - |
| Aq. HBr, sat. | В | А | А | С | - | А | А |
| Aq. HCl, sat. | А | А | А | В | А | А | А |
| Aq. ammonia solution | А | D | А | А | А | А | А |
| Aqueous caustic soda | А | В | А | А | В | А | А |
| Aqueous nitric acid | В | А | А | В | В | А | А |

*A: Very good resistance, B: Moderate resistance, C: poor resistance, D: very poor resistance

** A: Very good resistance, B: Moderate resistance, C: poor resistance


Please note: The resistance against the corresponding vapors is significantly better. Tabled values may vary by changing temperature and pressure.

| |
|------|
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |

| |
|------|
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |

For optimal distillation consider the BUCHI rule





For evaporation at 30 °C (boiling point)

- 1. instrument setting: 55 °C
- 2. cooling medium: 10 °C, 500 ml/min
- 3. recommended rotation setting: position 8-9
- 4. vacuum setting according to the Solvent List column 30 °C

Solvent List

| 6 | B | l | j | C | | i | | P | | |
|---|---|----|---|----|----|---|---|---|---|--|
| 5 | w | ÷Τ | z | ER | L. | A | N | D | / | |

| Solvent | Formula | Vacuum [mbar] for bp at: | | | |
|---|---|--------------------------|-------------|--|--|
| | 0110 | 30°C | 50°C | | |
| Acetic acid | $C_2H_4O_2$ | 26 | 72 | | |
| Acetone | C ₃ H ₆ O | 370 | 815 | | |
| Acetonitrile | C_2H_3N | 134 | 315 | | |
| n-Amylalcohol, n-Pentanol | C ₅ H ₁₂ O | 6 | 20 | | |
| <i>n</i> -Butanol | $C_4 H_{10}$ | 14 | 44 | | |
| <i>tert-</i> Butanol | $C_4H_{10}O$ | 78 | 231 | | |
| Chlorobenzene | C_6H_5CI | 22 | 56 | | |
| Chloroform | CHCl ₃ | 306 | 665 | | |
| Cyclohexane | $C_{6}H_{12}$ | 154 | 347 | | |
| 1,2-Dichloroethane | $C_2H_4CI_2$ | 137 | 315 | | |
| Dichloromethane | CH_2CI_2 | 685 | atm. press. | | |
| Diethylether | $C_4H_{10}O$ | 838 | atm. press. | | |
| trans-1,2-Dichloroethylene | C ₂ H ₂ Cl ₂ | 317 | 705 | | |
| Diisopropylether | $C_6H_{14}O$ | 251 | 545 | | |
| Dioxane | $C_4H_8O_2$ | 68 | 165 | | |
| Dimethylformamide (DMF) | C ₃ H ₇ NO | 6 | 17 | | |
| Ethanol | C_2H_6O | 97 | 276 | | |
| Ethylacetate | $C_4H_8O_2$ | 153 | 366 | | |
| Heptane | C_7H_{16} | 77 | 183 | | |
| Hexane | $C_{6}H_{14}$ | 241 | 525 | | |
| Isopropylalcohol | C ₃ H ₈ O | 78 | 231 | | |
| Isoamylalcohol | C ₅ H ₁₂ O | 9 | 29 | | |
| Methyl <i>tert-</i> buthyl ether (MTBE) | C ₅ H ₁₂ O | 413 | 835 | | |
| Methyl ethyl ketone (MEK) | C ₄ H ₈ O | 160 | 359 | | |
| Methanol | CH₄O | 236 | 607 | | |
| Pentane | $C_{5}H_{12}$ | 819 | atm. press. | | |
| <i>n</i> -Propanol | C ₃ H ₈ O | 37 | 115 | | |
| Pentachloroethane | C ₂ HCl ₅ | 8 | 21 | | |
| Tetrachloromethane | | 179 | 398 | | |
| Tetrahydrofurane (THF) | C_4H_8O | 234 | 539 | | |
| Toluene | C ₇ H ₈ | 48 | 118 | | |
| Trichloroethylene | C ₂ HCl ₃ | 119 | 275 | | |
| Water | H ₂ O | 42 | 120 | | |
| Xylene | $C_{8}^{2}H_{10}$ | 15 | 40 | | |

We are represented by more than 100 distribution partners worldwide. Find your local representative at

BÜCHI Labortechnik AG

CH – 9230 Flawil 1 T +41 71 394 63 63 F +41 71 394 65 65 buchi@buchi.com www.buchi.com

We are represented by more than 100 distribution partners worldwide. Find your local representative at

www.buchi.com

