


Protein is determined by the analysis of the nitrogen content. From this, the protein content is calculated. Protein consists of amino acids which contain nitrogen (N) in the amino group.



Other applications: TKN, ammonium, TVBN, urea, nicotine

Milling

The samples have to be homogeneous.

Micro	300 mL	500 mL
Reduction of chemicals by 80 %	Standard Kjeldahl	For high sample amounts

The actual weight of a sample depends on the nitrogen content as well as on the inhomogeneity of the sample.

Procedure:

- Select N % of sample
- Select titrant concentration
- Choose weight in order that the titrant consumption can be expected between 3 and 17 mL

- Conversion of K_2SO_4 to $KHSO_4$ (K_2SO_4 is a component of Kjeldahl Tablets)
approx. 2 – 3 mL
- Consumption by organic matter

Organic matter	H_2SO_4 g [mL]	Example: Salami	e.g. for 1.5 g weight (weight · organic matter)
Fat	9.7	27.3 %	$1.5 \cdot 9.7 \cdot \frac{27.3}{100} = 3.97$ mL
Protein	4.9	20.6 %	$1.5 \cdot 4.9 \cdot \frac{20.6}{100} = 1.51$ mL
Carbohydrates	4.0	0.0 %	$1.5 \cdot 4.0 \cdot \frac{0.0}{100} = 0.0$ mL
- Losses due to evaporation
approx. 1 mL/h
- Remaining volume
= Quantity of the used Kjeldahl Tablets (e.g. 10 g Tablets = 10 mL H_2SO_4)


 H_2SO_4 volume = conversion + (total consumption by org. matter) + evaporation + remaining volume
 3 mL + (3.97 + 1.51 + 0.00) mL + 1 mL + 10 mL = 18.48 mL ~18 mL

Article	Order code	Composition	Weight
Titanium	11057980	3.5 g K_2SO_4 / 0.105 g $CuSO_4 \cdot 5 H_2O$ 0.105 g TiO_2	3.71 g
Missouri	11057982	4.98 g K_2SO_4 / 0.02 g $CuSO_4 \cdot 5 H_2O$	5 g
ECO	11057983	3.998 g K_2SO_4 / 0.002 g $CuSO_4$	4 g
Antifoam	11057984	0.97 g Na_2SO_4 / 0.03 g Silicon Antifoam	1 g
Titanium Micro	11057981	1.5 g K_2SO_4 / 0.045 g $CuSO_4 \cdot 5 H_2O$, 0.045 g TiO_2	1.59 g
Copper Micro	11057985	1.5 g K_2SO_4 / 0.15 $CuSO_4 \cdot 5 H_2O$	1.65 g

The aim of Kjeldahl Tablets is to accelerate the digestion by:

- Catalysis by metal salts
- Increasing the boiling point of H_2SO_4 by sulfate salts (K_2SO_4)

Attention: With boiling temperatures above 390 °C nitrogen losses are possible \Rightarrow leads to too low results!


$$\text{Digestion: } (C_n H_n N_n O_n) + H_2SO_4 \rightarrow CO_2 + SO_2 + H_2O + NH_4^+$$


Organic matter is destroyed by boiling in concentrated sulfuric acid. Kjeldahl Tablets raise the boiling point and accelerate the process.




Distillation: $\text{NH}_4^+ + \text{OH}^- \rightleftharpoons \text{NH}_3 (\text{Gas}) + \text{H}_2\text{O}$

The digestion mixture is alkalinized with NaOH prior to distillation to free up the ammonia. The ammonia is steam distilled into an acidic receiver solution.

The titration



The pH in the acidic receiver solution rises upon addition of ammonia. The nitrogen and protein content is then determined by titration of the borate complex.



Endpoint titration boric acid: $\text{B(OH)}_3 + 2 \text{H}_2\text{O} \rightleftharpoons \text{B(OH)}_4^- + \text{H}_3\text{O}^+$ ($\text{pK}_a = 9.24$)

The result calculation

$$\%N = \frac{[V(1) - V(B)] \cdot F \cdot c \cdot f \cdot M(N)}{m \cdot 1000} \cdot 100 \%$$

% P = % N · PF

V(1):	consumption of titrant, sample [mL]
V(B):	average consumption of titrant, blank [mL]
F:	molar reaction factor (1 = HCl, 2 = H ₂ SO ₄)
c:	concentration of titrant [mol / L]
f:	factor of titrant
M(N):	molecular weight of N (14,007 [g/mol])
m:	sample weight [g]
1000:	conversion factor (mL in L)
PF:	protein factor
% N:	% of weight of N
% P:	% of weight of protein

Example:		
Titrant:	0.25 mol/L H_2SO_4	$(7.5 \text{ mL} - 1 \text{ mL}) \cdot 2 \cdot 0.25 \cdot 1 \cdot 14.007 \cdot 100\%$
	$f = 1.000$	$\frac{}{0.750 \text{ g} \cdot 1000} = 6.07\% \text{ N}$
Weighing:	0.750 g	
Blank value:	1.0 mL	
Titration volume for sample: 7.5 mL		$6.07\% \text{ N} \cdot 6.25 = 37.94\% \text{ P}$

Food	Amount N(%) protein factor		Food	Amount N(%) protein factor	
Egg / meat / fish	16.0	6.25	Almonds	19.3	5.18
Milk / ice cream / casein	15.7	6.38	Peanuts / brasil nuts	18.3	5.46
Wheat / baked prod. / oat	17.5	5.70	Tee nut / coconut	18.9	5.30
Rye / soy	17.5	5.70	Rice	19.0	5.26
Barley / maize / pulses	16.0	6.25	Beer / brewing sugars	16.0	6.25