


FatExtractor E-500

Technical data sheet

The declaration of the total fat or crude fat is required for most foods and feeds. The FatExtractor E-500 is designed for solid-liquid extraction for quick and compliant fat determination. The design of the glass assemblies and the high-speed heaters combined with sophisticated process control allow for the fastest and most reproducible extraction processes with full compliance.

By simply changing the glass assembly, the FatExtractor E-500 complies with standard methods such as Soxhlet, Hot Extraction (HE) or Twisselmann (ECE). The glass assembly is interchangeable for a later conversion to another extraction method. Suitable for extractions with petroleum ether, diethyl ether, hexane or chloroform.





Description of function

The FatExtractor E-500 is designed to carry out the following solid-liquid extraction methods:

- Soxhlet Extraction (see Chapter "Soxhlet Extraction")
- Hot Extraction (see Chapter "Hot Extraction")
- Economic Continuous Extraction (see Chapter "Economic Continuous Extraction")

Soxhlet Extraction

Step 1: Extraction

- The sample is located in the extraction chamber.
- The beaker contains the solvent.
- The solvent is heated, vapor rises up to the condenser, condenses and drops into the extraction chamber with the sample.
- The magnetic valve is closed, the solvent is collected up to the optical sensor and extracts the analyte.
- When the optical sensor is reached, the magnetic valve opens and the solvent containing the analyte flows back into the beaker.

Step 2: Rinsing

- The solvent is heated, vapor rises up to the condenser, condenses and drops into the extraction chamber with the sample.
- The magnetic valve is open, the solvent flows back into beaker, the solvent is not collected.

Step 3: Drying

- The solvent is heated, vapor rises up to the condenser, condenses and flows into tank.
- The analyte remains in the beaker.

Hot Extraction

Step 1: Extraction


- The sample is located in the beaker.
- The beaker contains the solvent.
- The solvent is heated, vapor rises up to the condenser, condenses and drops into the beaker with the sample.


Step 2: Rinsing

- The solvent in the beaker is heated up and evaporated.
- The vapor rises up to the condenser.
- The condensed solvent flows into the beaker with the sample.
- The tank bottle valve opens periodically and condensed solvent flows in the tank bottle.
- The solvent level decreases.

Step 3: Drying

- The solvent is heated, vapor rises up to the condenser, condenses and flows into tank.
- The analyte remains in the beaker.





Economic Continuous Extraction

Step 1: Extraction

- The sample is located in the extraction chamber.
- The beaker contains the solvent.
- The solvent is heated, vapor rises up around the sample to the condenser, condenses and drops back into extraction chamber through the sample into the beaker.

Step 2: Drying

- The solvent is heated, vapor rises up around the sample to the condenser, condenses and flows into tank.
- The analyte remains in the beaker.




Order code FatExtractor E-500

The FatExtractor E-500 can be ordered in the following configurations: Soxhlet, Hot Extraction (e.g. Randall) or Economic Continuous Extraction (e.g. Twisselmann). For large sample volumes, use the Soxhlet LSV version, suitable for extraction thimbles with a diameter of 43 mm.

The FatExtractor E-500 can be upgraded with the Analyte protection sensor for increased user safety and to prevent any loss of analyte during the extraction process. The Pro Interface consists of a 7 inch color touch display for even more intuitive use and overview on the interface.

For environmentally friendly operation the FatExtractor E-500 can be bundled with Recirculating Chiller F-305/F-308 for saving cooling water. The FatExtactor E-500 has to be operated either with a recirculating chiller or a tap water valve.

Choose the configuration according to your needs:




Order code FatExtractor E-500 System

The FatExtractor E-500 can be bundled with the HydrolEx H-506 for total fat determination and a Recirculating Chiller F-305/F-308.

HydrolEx H-506 performs convenient and safe acid hydrolysis as a step prior to extraction. It is an essential work step of the total fat determination where matrix structures enclosing the fat fraction of food and feed samples are broken up.

For environmentally friendly operation the FatExtractor E-500 can be bundled with Recirculating Chiller F-305/F-308 for saving cooling water. The FatExtactor E-500 has to be operated either with a recirculating chiller or a tap water valve.

Choose the configuration according to your needs:



Scope of delivery
All configurations are supplied ready to use.

	Hot Extraction	Economic Continuous Extraction	Soxhlet Extraction	Soxhlet Extraction LSV
FatExtractor E-500	1	1	1	1
Condenser E-500	6	6	6	6
Extraction glass chamber ECE	-	6	-	-
Extraction glass chamber Soxhlet	-	-	6	-
Extraction glass chamber Soxhlet LSV	-	-	-	6
Soxhlet assembly cpl.	-	-	6	6
Extraction beaker	-	6	6	6
Extraction beaker HE	6	-	-	-
Sealing E-500	6	12	12	12
Set of glass sample tube holder (3 pcs.)	2	2	2	-
Set of glass sample tube holder LSV (3 pcs.)	-	-	-	2
Set of holders for thimbles 25 mm (3 pcs.)	2	2	2	-
Set of holders for thimbles 33 mm (3 pcs.)	2	2	2	2
Set of holders for thimbles 43 mm (3 pcs.)	-	-	-	2
Extraction thimbles 25 x 100 mm	6	6	6	-
Extraction thimbles 33 x 94 mm	6	6	6	6
Extraction thimbles 43 x 118 mm	-	-	-	6
Solvent tank cpl.	1	1	1	1
Cooling water hose 3 m	2	2	2	2
Beaker tong	1	1	1	1
Extraction beaker carrier	-	1	1	1
Extraction beaker carrier HE	1	-	-	-
Funnel	1	1	1	1
Power cord	1	1	1	1

	Hot Extraction	Economic Continuous Extraction	Soxhlet Extraction	Soxhlet Extraction LSV
Pliers for glass sam- ple tube with frit	1	1	1	1
Operation manual	1	1	1	1

Technical data

FatExtractor E-500

Specification	E-500
Power consumption	1300 W
Connection voltage	100 - 240 ± 10 % VAC
Fuse	10 A
Frequency	50 / 60 Hz
Overvoltage category	II
Pollution degree	2
Protection Class	1
Dimensions (W x D x H) (with glassware Soxhlet Extraction)	638 x 595 x 742 mm
Dimensions (W x D x H) (with glassware Hot Extraction)	638 x 595 x 613 mm
Dimensions (W x D x H) (with glassware Economic Continuous Extraction)	638 x 595 x 622 mm
Minimum clearance (W x D)	200 mm
Minimum clearance (H)	500 mm
Weight (without glassware)	41.8 kg
Weight (with glassware Soxhlet Extraction)	49.6 kg
Total heating power (rated)	600 W
Total heating power (maximum)	1200 W
Hose connection	6 / 9 mm
Minimum water flow	100 mL/min
Allowed water pressure (nominal value)	6 bar
Allowed water pressure (maximum)	8 bar
Inlet cooling medium temperature	25 °C below the boiling point of the solvent
Number of extraction positions	6
Solvent tank volume	2 L

Specification	E-500
Max. filling level	120 mL
(Extraction glass chamber Soxhlet)	
Max. filling level	195 mL
(Extraction glass chamber Soxhlet LSV)	
Max. working volumes	175 mL
(Beaker)	
Max. working volumes	100 mL
(Beaker HE)	
Language	DE, EN, IT, ES, FR, JA, CN, PL,
	RU
Method storage	20 methods
(Interface)	
Method storage	40 methods
(Pro interface)	

Ambient conditions

For indoor use only.

Max. altitude above sea level	2000 m
Ambient temperature	5–40 °C
Maximum relative humidity	80% for temperatures up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C
Storage temperature	max. 45 °C

Conversion and upgrade kits

The FatExtractor E-500 can be converted into another configuration with an easy change of glass assembly. Choose the needed conversion kit according to the matrix.

from	Soxhlet Extraction	Hot Extraction	Economic Continuous Extraction	Soxhlet Extraction LSV
Soxhlet Extraction		11068489	11068492	
Hot Extraction	11068487		11068487	11068487
Economic Continuous Extraction	11068488	11068491		11068488
Soxhlet Extraction LSV		11068490	11068493	

	Order no.
Conversion kit from Soxhlet / ECE to HE	11068487
incl.: 6 beakers HE (3 \times 11067475) , beaker brackets HE (2 \times 11067829), set of draining tubes for HE (11067480), beaker carrier cover for HE	
Conversion kit from Soxhlet / Universal to ECE	11068488
Includes 6 Extraction glass chamber ECE (11062499), set of draining tubes for ECE (11067479)	
Conversion kit from HE to Soxhlet	11068489
Includes 6 Extraction glass chamber Soxhlet (11062496), 6 beaker (3 x 11067474), 6 beaker brackets (2 x 11067828), 6 Soxhlet assemblies cpl. (11067065) and a chamber rack (11067077), set of draining tubes for Soxhlet (11067478), two safety shields top, cpl. (2 x 11067832)	
Conversion kit from HE to ECE	11068491
Includes 6 Extraction glass chamber ECE (11062499), 6 beaker (3 x 11067474), 6 beaker brackets (2 x 11067828), set of draining tubes for ECE (11067479) and a chamber rack (11067077)	
Conversion kit from ECE to Soxhlet	11068492
Includes 6 Extraction glass chamber Soxhlet (11062496), 6 Soxhlet assemblies cpl. (11067065) and set of draining tubes for Soxhlet (11067478), two safety shields top, cpl. (2 x 11067832)	
Conversion kit HE to Soxhlet LSV	11068490
Includes 6 Extraction glass chamber Soxhlet LSV (11062497), 6 beaker (3 \times 11067474), 6 beaker brackets (2 \times 11067828), 6 Soxhlet assemblies cpl. (11067065) and a chamber rack (11067077), set of draining tubes for Soxhlet (11067478), two safety shields top, cpl. (2 \times 11067832)	
Conversion kit ECE to Soxhlet LSV	11068493
Includes 6 Extraction glass chamber Soxhlet LSV (11062497), 6 Soxhlet assemblies cpl. (11067065) and a set of draining tubes for Soxhlet (11067478), two safety shields top, cpl. (2 x 11067832)	

	Order no.
Upgrade kit Analyte protection	11068524
For later upgrade of a FatExtractor E-500 with analyte protection sensor	
Upgrade kit Interface Pro	11068525
For later upgrade of a FatExtractor E-500 with a Pro interface (7 " touch screen)	

Chiller

No. Units	Ambient Temperature	Chiller
1	< 30 °C	Recirculating Chiller F-305
1	< 40 °C	Recirculating Chiller F-308
2	< 30 °C	Recirculating Chiller F-308
2	< 40 °C	Recirculating Chiller F-314

	Order no.
Recirculating Chiller F-305	11F30501
550 W at 15 °C, Display, 230 V	
Recirculating Chiller F-305	11F30502
550 W at 15 °C, Display, 115 V	
Recirculating Chiller F-308	11F30801
900 W at 15 °C, Display, 230 V	
Cooling capacity 900 W at 15 °C, for temperatures from -10 to 25 °C	
Recirculating Chiller F-308	11F30802
900 W at 15 °C, Display, 115 V	
Cooling capacity 900 W at 15 °C, for temperatures from -10 to 25 °C	
Recirculating Chiller F-314	11F31401
1400 W at 15 °C, Display, 230 V	
Cooling capacity 1400 W at 15 °C, for temperatures from -10 to 25 °C	
Recirculating Chiller F-314	11F31402
1400 W at 15 °C, Display, 115 V	
Cooling capacity 1400 W at 15 °C, for temperatures from -10 to 25 °C	

Spare parts

	Order no.	Image
Extraction glass chamber Soxhlet	11062496	
Extraction glass chamber Soxhlet LSV	11062497	

	Order no.	Image
Extraction glass chamber ECE	11062499	
Set of beakers, 2 pcs.	11067474	
Set of beakers HE (2pcs)	11067475	
Condenser E-500 cpl.	11067063	
Condenser flange E-500	11067817	
Condenser tank bottle	11065966	
Tank adapter, PTFE	11064590	
Tank bottle 2 L, GL 45	11070509	
Set of gliding elements including magnets, 10 pcs.	11067827	

	Order no.	Image
Soxhlet assembly cpl.	11067065	
One part constiting of magnetic valve and level sensor for extraction glass chamber Soxhlet		
Protection shield top, cpl.	11067832	
Protection shield bottom, cpl.	11067831	
Set of seals E-500, FKM, 6 pcs.	11069012	
Membrane with anchor for magnetic valve unit	037534	
Joint clip	11070136	
Set of beaker brackets, 3 pcs.	11067828	
Set of beaker brackets HE, 3 pcs.	11067829	
Reflectorfoil analyte protection, 6 pcs.	11068522	
Knurled nut and sealing disks, 2 pcs. Connection of the draining tubes to the condensers	11068520	
Set of flange locks, 6 pcs.	11067833	
Silicone hose D6/9 L=3 m	048355	
Set of draining tubes SOX, FEP	11067478	
The draining tubes connect the receiving funnel in the condensers with the tank valve to drain solvent into the tank.		
Set of draining tubes HE, FEP	11067480	
The draining tubes connect the receiving funnel in the condensers with the tank valve to drain solvent into the tank.		
Set of draining tubes ECE, FEP	11067479	

Accessories

	Order no.	Image
Holder for glass sample tubes, stainless steel	11067219	

	Order no.	Image
Holder for glass sample tubes, PTFE	11067220	
Holder for extraction thimbles (diameter 25 - 43 mm)	11068443	
Extraction beaker carrier	11067042	
Allows to carry 6 beakers (11067474)		
Extraction beaker HE carrier Allows to carry 6 beakers HE (11067475)	11067493	
Set condenser insulations E-500, 6 pcs.	11069078	
The insulation of the condensers prevent condensing water and is recommended in high humidity environment		
Set insulation cooling water hoses	11069079	
The insulation of the water hoses prevent condensing water and is recommended in high humidity environment.		
Support solvent supply	11068306	
Allows to fix the tubes of solvent dispensers to the condensers for convenient solvent addition.		
Cooling water valve. 24VAC	031356	
Valve opens cooling water feed during distillation. Meant to be used with a vacuum controller/interface.		
Extension cable for the cooling water valve, 2 m	11069477	
Turning disk	11067985	
Allows for turning the instrument for easier access.		
StatusLight cpl., incl. communication cable	11068959	
Indicates the status of the instrument (instrument is ready to use, has an error or is in operation).		

	Order no.	Image
Funnel	11067473	

Consumables

	Order no.
Quartz sand 0.3 - 0.9 mm, 2.5 kg	037689
Celite® 545, 1 kg	11068920

Glass sample tubes and extraction thimbles

	Order no.	Image
Glass sample tubes with frit, 6 pcs.	11067497	
Working volume: 64 mL		
Filling volume: 82 mL		
Glass sample tubes LSV with frit, 6 pcs.	11067814	
Working volume: 116 mL		
Filling volume: 144 mL		
Extraction thimbles 25 x 100 mm, 25 pcs.	018105	
Working volume: 44 mL		
Extraction thimbles 33 x 94 mm, 25 pcs.	11058983	
Working volume: 64 mL		
Extraction thimbles, Set. 25 pcs, 43 x 118 mm, cellulose	018106	
For Soxhlet extraction unit.		
Working volume: 150 mL		

Holder for extraction thimbles

	Order no.
Set of holders for glass sample tubes with frit, PTFE, 3 pcs.	11067485
Set of holders for LSV glass sample tubes, PTFE, 3 pcs.	11067486
Holders for thimbles d25, PTFE, 3 pcs.	11067488
Holders for thimbles d33, PTFE, 3 pcs.	11067490
Holders for thimbles d43, PTFE, 3 pcs.	11067491

Documentation

	Order no.
IQ/OQ for FatExtractor E-500	11068557
Installation/Operation qualification for the FatExtractor E-500 incl. binder with documents IQ/OQ, inspection and basic tag (English)	
Repeating OQ for FatExtractor E-500	11068559
Repeating operation qualification for FatExtractor E-500, incl. document, inspection and basic tag (English)	