



Ultrasonic Package

Maximize your particle size


Generate larger particles in laboratory scale spray drying. The Ultrasonic Package allows an easier handling of your product and increases its flowability.



Improved powder handling
Improved dosing and handling due to better flow ability



Highly performant
Large particles with narrow size distribution (10 – 60 μm)



Easy operation
Easy installation and operation of the Ultrasonic Package and the Maxi Spray Dryer B-290

Ultrasonic Package: Your most important benefits



improved powder handling

- Large particles (10 – 60 μm)
- Better flowability of the powder for Improved handling
- Easy scale up process due to particle size close to industrial scale
- Better distribution of active materials within powder



Highly performant

- Narrow particle shape and size distribution
- Saving valuable material (> 5 mL) thanks to small sample amount
- Reproducible powder production
- High yields (up to 70 %)
- Shorter time to optimize formulations



Easy operation

- Easy upgrade of the Mini Spray Dryer B-290
- Intuitive handling
- Fast setup and short cleaning times
- Visible process due to glass assembly

Complete your portfolio



**Mini Spray Dryer
B-290**

World leading laboratory Spray Dryer



**Nano Spray Dryer
B-90**

Spray Dryer for small samples and particles



**Inert Loop
B-295**

Spray Drying organic solvents



**Dehumidifier
B-296**

Reproducible drying air for Spray Drying

